



Monocular Student Microscope A11.1506-A1 Biological Vertical Head Tube

Our Product Introduction

for more products please visit us on cnoec.com

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A11.1506
- Minimum Order Quantity: 1 pc
- Price: Negotiation
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Total Magnification: 40x-400x
- Head: Monocular Vertical Head Tube
- Eyepiece: WF10x
- Objectives: 4x, 10x, 40xs
- Focusing: Coarse Adjustment Range: 20mm
- Stage: Plain Stage With Clips, 100mmx90mm
- Illuminator: Concave-Plano Mirror Dia. 40mm
- Highlight: **binocular biological microscope, biological compound microscope**



More Images



Product Description

OPTO-EDU A11.1506-A1 Monocular Cheapest professional educational student microscope

Total Magnification: 40X~400X

Monocular Vertical Tube, Coarse Focus

Simple Structure And Nice Appearance

Featured By Clear Image, Easy Operation

The Biological Microscope Model A11.1506, Can Be Used In Teaching Experiment For Primary And Secondary Schools. It Is Featured By Clear Image, Easy Operation, Simple But Reliable Structure And Nice Appearance. Stration And Other Occasions.

A11.1506 Specification		A1	A2	A3
Total Magnification	40x-400x	●	●	●
Head	Monocular Vertical Head Tube	●	●	-
	Monocular 45°Inclined,360°Rotatable	-	-	●
Eyepiece	WF10x	●	●	●
Objectives	4x	●	●	●
	10x	●	●	●
	40xs	●	●	●
Focusing	Coarse Adjustment Range: 20mm	●	-	-
	Seperate Coarse And Fine Adjustment,Coarse Adjustment Range: 20mm, Fine Adjustment Range : 1.4~1.6mm	-	●	●
Stage	Plain Stage With Clips,100mmx90mm	●	●	●
Illuminator	Concave-Plano Mirror Dia. 40mm	●	●	●
Optional Accessories				
Illuminator	Electrical Illuminator (110V/220V 20W)			

OPTO-EDU

Biological Microscope

A11.1506



 Total Magnification: 40X~400X

 Monocular Vertical Tube, Coarse Focus

 Simple Structure And Nice Appearance

 Featured By clear Image, Easy Operation




Don't Miss Best Microscope For You!

A11.1506 BIOLOGICAL MICROSCOPE



 A11.1506-A1 Monocular Coarse Focus Biological Microscope




 A11.1506-A2 Monocular Separate Coarse And Fine Focus Adjustment Biological Microscope



▲ A11.1506-A3 Inclined Monocular Separate Coarse And Fine Focus Adjustment Biological Microscope

USER INSTRUCTION



01 WORKING STAGE

Put The Sample Slice Required To Observe On The Center Of The Stage, At The Center Of The Round Hole. Press Firmly With Two Pieces Of Spring.

02 ILLUMINATOR

Face The Mirror Toward A Light Source, Such As A Fluorescent Lamp Or A Window, To Conduct Light Into The Microscope And Obtain A Bright Field Of View.



03 OBJECTIVE

Achromatic 4x, 10x, 40x Objectives Included. Put The Objective Needed To Observe To The Center Under Low Magnification 4x, Rotate Coarse Focus Till Presence Of The Clearest Image Before High Magnification Observation.

04 FOCUSING

After The Replacement Of The High Magnification Lens, View From Eyepiece, In Case Of An Unclear Image. Rotate Coarse Focus Until A Clear Image Is In The View.



PRODUCT DESCRIPTION



PACKING AND SHIPPING

Packing : Cardbox And Polyfoam + Strong Master Carton 15*21*34cm, G.W.: 2.0Kgs

Shipping Time : 3-10 Days From Stock



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China