



CCD 40x - 1000x Portable Digital Optical Microscope A31.1006 For Laboratory

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE,
- Model Number: A31.1006
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification


- Applications: Laboratory
- Certification: CE
- Drawtube: Monocular/Binocular
- Illumination: LED Brightness Adjustable
- Max Magnification: 1000X
- Theory: Digital Optical Microscope
- Highlight: **digital usb microscope,**
optical microscope with digital camera

Product Description

Monocular/Binocular CCD 40x - 1000x Portable Digital Microscope A31.1006 For Laboratory

A31.1006 Digital Biological Microscope		A	B
Head	Monocular head, Inclined at 45°, 360°rotatable	●	
	Binocular head, Inclined at 45°, 360°rotatable		●
Eyepiece	WF10x/18		
Objective	Achromatic 4x,		
	Achromatic 10x		
	Achromatic 40x		
	Achromatic 100x		
Nosepiece	Quadruple Nosepiece		
Stage	Double Layers Mechanical Stage 115×125mm/70×20mm		
Condenser	Abbe N.A.=1.25 with iris diaphragm(Φ2-Φ28mm)and frosted filter		
Focusing	Coaxial coarse and fine adjustment		
Camera System	Valid Pixel 1280×1024(1.3M Pixel)		
	Output Mode USB2.0		
	Operation System WINDOWS 2000/XP/VISTA 32 Bit or 64 Bit		
	Software ScopelImage 9.0		
Illumination	1W Powerful LED Illumination, brightness adjustable		





Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China