



A22.1504 Binocular Stereo Zoom Microscopes 20x-40x Magnification With 57mm Working Distance

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A22.1504
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Magnification: 20x-40x
- Applications: Medical|Clinic
- Certification: CE|Rohs
- Drawtube: Binocular
- Product Name: Stereo Microscope
- Working Distance: 57mm
- Highlight: **stereo zoom microscopes,
stereo dissecting microscope**

Product Description

Binocular Stereo Microscope 20x-40x Magnification A22.1504 With 57mm Working Distance

Specifications:

A22.1504 Specification								
Model	Eyepiece	Objective	Total Magnification	Coarse Adjustment Range	Working Distance	Plate	Illuminator	Light Source
3A	WF10X	2x/4x	20x-40x	40mm	57mm	White/Black	No	No
3B	WF10X	2x/4x	20x-40x	40mm	57mm	White/Black	Incident	Halogen Lamp 12V/10W
3C	WF10X	2x/4x	20x-40x	40mm	57mm	White/Black	Incident & Transmitted	Halogen Lamp 12V/10W



A22.1204-1A

www.cnoec.com



A22.1204-1B



A22.1204-1C



Opto-Edu (Beijing) Co., Ltd.



0086 13911110627



sale@optoedu.com



cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China