



Led Fluorescent Light Microscopy A16.2601-I With Double Layer Mechanical Stage

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A16.2601-L
- Minimum Order Quantity: 1 pc
- Price: Negotiation
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Eyepiece: WF10x/20mm
- Focusing System: CE|Rohs
- Nosepiece: Quintuple
- Stage: Double Layer Mechanical, Size 175x145mm, Moving Range 78x55mm, Two Slide Holder
- Condenser: Abbe N.A.1.25 With Aperture Iris Diaphragm
- Transmitted Illumination: Koehler, 6V20W Halogen Lamp With Brightness Adjustable
- Highlight: **inverted fluorescence microscope, epifluorescence microscope**

Product Description

LED Fluorescence Light Microscope A16.2601-L

Specification	A16.2601-L-B2	A16.2601-L-T2	A16.2601-L-B4	A16.2601-L-T4
	Binocular	Trinocular	Binocular	Trinocular
Seidentopf Viewing Head	30° inclined, Interpupillary distance 50-75mm, For trinocular head, light distribution (binocular tube:photo tube): both 100:0 and 20:80		30° inclined, Interpupillary distance 50-75mm, For trinocular head, light distribution (binocular tube:photo tube): both 100:0 and 20:80	
Eyepiece	WF10x/20mm			
Objective	Infinity Plan Achromatic Objective 4x / 0.10			
	Infinity Plan Achromatic Objective 10x / 0.25			
	Infinity Plan Achromatic Objective 20x / 0.40(s)			
	Infinity Plan Achromatic Objective 40x / 0.66(s)			
Focusing System	Infinity Plan Achromatic Objective 100x / 1.25(s, oil)			
Nosepiece	Low positioned, coaxial coarse and fine focusing knobs			
Stage	Quintuple			
Condenser	Double layer mechanical, size 175x145mm, moving range 78x55mm, two slide holder			
	Abbe N.A.1.25 with aperture iris diaphragm			
Fluorescent Attachment	Fluorescent Filter (B,G)		Fluorescent Filter (UV,V) are under development	
	Six-hole Disc Epi-Fluorescence Set			
	3W LED Lamp House			
	3W LED Power Supply			
Transmitted Illumination	Koehler, 6V20W Halogen Lamp With Brightness Adjustable			
Filter	Blue			
Power Supply	100-240V, 50/60Hz			
Packing	Carton, Size 99x34x40cm, G/N weight 23/20Kgs			
Optional Accessories			Item No.	
Eyepiece	WF16x/13mm		A51.2602-16	
	WF10x/18mm with pointer		A51.2604-10	
	WF10x/20mm with micrometer 0.1mm/Div		A51.2605-10A	
Objective	Infinity Plan Objectives 2.5x/0.07		A52.2606-2.5	
	Infinity Plan Objectives 60x/0.80(s)		A52.2606-60	
Light Source	Upgrade to 6V 30W		A56.2650	
LED				
Fluorescent Attachment	Chroma Filter, Wavelength 455nm Speically Used for Tuberculosis Inspection		A5F.2602-TB	
Condenser	Swing out achromatic condenser N.A.0.9/0.13		A56.2614-2	
Filter	green,		A56.2616-G	
	gray		A56.2616-N	
	yellow		A56.2616-Y	
CCD Adapter	C-Mount 0.57x		A55.2601-BK057	
	C-Mount 1.0x		A55.2601-BK10	
Digital Camera Adapter	For CANON digital camera		A55.2621	
	For OLYMPUS digital camera		A55.2620	
	For NIKON SLR digital camera		A55.2622-N	
	For SONY SLR digital camera		A55.2622-S	
	For Canon SLR digital camera		A55.2622-C	



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China