



## 0.7 - 4.5x Pole Stand Zoom Stereo Optical Microscope Opto-edu Binocular A23.3645-B1

Our Product Introduction

for more products please visit us on [cnoec.com](http://cnoec.com)

### Basic Information

- Place of Origin: China
- Brand Name: OPTO-EDU
- Certification: CE, Rohs
- Model Number: A23.3645-B1
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



### Product Specification

- Eyepiece: WF10x/20mm
- Zoom Lens: 0.7~4.5x
- Zoom Ratio: 1:6.5
- Magnification: 7x-45x, Up To 3.5x-180x With Optional Eyepieces And Auxiliary Objectives
- Working Distance: 100mm, With Auxiliary Objectives 30~165mm
- Field Of View: 1~95.2mm
- Stand: Pole Stand, No Light Source
- Highlight: **stereo zoom microscopes,  
stereo dissecting microscope**



### More Images



## Product Description

0.7-4.5x, Binocular, WF10x/20, Magnification 7~45x, Optional Eyepieces And Auxiliary Objectives Can Up To 3.5-270x, B1 Pole Stand + A1 Holder, Stand Size 200\*255\*22mm, Pole Height 240mm, Dia.32mm, No Light Source, Pressure plate and Dia.95mm black and white board



## PRODUCT DESCRIPTION

A23.3645-B1 Zoom Stereo Microscope Specification		B	T
Head	Binocular Head, 45° Inclined, 360°Rotatable, ±5 Diopter Adjustable, Interpupillary Adjustment 54-75mm	O	
	Trinocular Head, 45° Inclined, 360°Rotatable, ±5 Diopter Adjustable, Interpupillary Adjustment 54-75mm, Eyepiece Tube Adapter		O
Eyepiece	WF10x/20mm		
Zoom Lens	0.7~4.5x		
Zoom Ratio	1:6.5		
Magnification	7x-45x, Up To 3.5x-270x With Optional Eyepieces And Auxiliary Objectives		
Working Distance	100mm, With Auxiliary Objectives 30~165mm		
Field Of View	1~95.2mm		
Stand	Pole Stand, No Light Source		
Optional Accessories		Cata.No.	
Head	Binocular Head + WF10x/20mm	A53.3645-B	
	Trinocular Head + WF10x/20mm	A53.3645-T	
Eyepiece	WF10x/20mm, High-eyepoint	A51.3621-1020	
	WF10x/22mm, High-eyepoint	A51.3621-1022	
	WF10x/23mm, High-eyepoint	A51.3621-1023	
	WF15x/15mm, High-eyepoint	A51.3621-1515	
	WF20x/10mm, High-eyepoint	A51.3621-2010	
	WF25x/9mm, High-eyepoint	A51.3621-2509	
	WF30x/8mm, High-eyepoint	A51.3621-3008	
	WF10x/22mm, High-eyepoint, Diopter Adjustable	A51.3622-1022	
	WF15x/15mm, High-eyepoint, Diopter Adjustable	A51.3622-1515	
	WF20x/10mm, High-eyepoint, Diopter Adjustable	A51.3622-2010	
	WF25x/9mm, High-eyepoint, Diopter Adjustable	A51.3622-2509	
	WF30x/8mm, High-eyepoint, Diopter Adjustable	A51.3622-3008	
Auxiliary Lens	0.5X,W.D. 165mm	A52.3621-0.5x	
	1.5X,W.D. 45mm	A52.3621-1.5x	
	2.0X,W.D. 30mm	A52.3621-2.0x	

CCD Adapter	1.0x, For CTV	A55.3621-1.0x
	0.50x, For 1/2" CCD	A55.3621-0.5x
	0.35x, For 1/3" CCD	A55.3621-0.35x
	0.50x, For 1/2" CCD	A55.3621-N0.5x

## PRODUCT DISPLAY



● Eyepiece WF10X/20mm

● 0.7x-4.5x, Magnification Up To 3.5x-180x With Optional Eyepieces And Auxiliary Objectives




A54.3640-A1



A54.3630-B1



● B1 Pole Stand + A1 Holder




A23.3645 Zoom Stereo Microscope Optical Data Sheet								
Objective	Standard 1.0x		Auxiliary 0.5X		Auxiliary 1.5X		Auxiliary 2X	
W.D.	100mm		165mm		45mm		30mm	
Eyepiece	Mag	View Field	Mag	View Field	Mag	View Field	Mag	View Field
10X/20mm	7.0X	28.6mm	3.5X	57.1mm	10.5X	19mm	14X	14.3mm
	45.0X	4.4mm	22.5X	8.9mm	67.5X	3mm	90X	2.2mm
15X/15mm	10.5X	21.1mm	5.25X	42.8mm	15.75X	14.3mm	21X	10.7mm
	67.5X	3.3mm	33.75X	6.7mm	101.25X	2.2mm	135X	1.7mm
20X/10mm	14.0X	14.3mm	7.0X	28.6mm	21.0X	9.5mm	28X	7.1mm
	90.0X	2.2mm	45.0X	4.4mm	135.0X	1.5mm	180X	1.1mm

## PACKING AND SHIPPING

**Packing** : By Polyfoam + Strong Master Carton 39\*29\*41cm, G.W.: 5.5Kgs

**Shipping Time** : 3-10 Days From Stock



**Opto-Edu (Beijing) Co., Ltd.**



0086 13911110627



sale@optoedu.com



cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China