



Trinocular Stereo Optical Microscope With Optional Eyepieces / Auxiliary Objectives

Our Product Introduction

for more products please visit us on cnoec.com

Basic Information

- Place of Origin: China
- Brand Name: OPTO-EDU
- Certification: CE, Rohs
- Model Number: A23.3645N-B1T
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Eyepiece: WF10x/20mm
- Zoom Lens: 0.7~4.5x
- Zoom Ratio: 1:6.5
- Magnification: 7x-45x, Up To 3.5x-180x With Optional Eyepieces And Auxiliary Objectives
- Working Distance: 100mm, With Auxiliary Objectives 30~165mm
- Field Of View: 1~95.2mm
- Stand: Pole Stand, No Light Source
- Highlight: **stereo zoom microscopes,**
stereo dissecting microscope



Product Description

Hinged triad, 50/50%, 20/80%, WF10X/20mm High Eyepieces, With A Variety Of Base Options, Up to 3.5-270x With Optional Eyepieces And Auxiliary Objectives.

A23.3645N-B1 Zoom Stereo Microscope Specification		B	T
Head	Binocular Head, 45° Inclined, 360°Rotatable, ±5 Diopter Adjustable, Interpupillary Adjustment 54-75mm	O	
	Trinocular Head, 45° Inclined, 360°Rotatable, ±5 Diopter Adjustable, Interpupillary Adjustment 54-75mm, Eyepiece Tube Adapter		O
Eyepiece	WF10x/20mm, Both With Eye-cup		
Zoom Lens	0.7~4.5x		
Zoom Ratio	1:6.5		
Magnification	7x-45x, Up To 3.5x-180x With Optional Eyepieces And Auxiliary Lens		
Working Distance	100mm, With Standard 1x Lens		
Field Of View	1~95.2mm		
Stand	Pole Stand, No Light Source		
Optional Accessories		Cata.No.	
Head	Binocular Head + WF10x/20mm	A53.3645N-B	
	Trinocular Head + WF10x/20mm	A53.3645N-T	
Eyepiece	WF10x/20mm, High-eyepoint	A51.3621-1020	
	WF10x/22mm, High-eyepoint	A51.3621-1022	
	WF10x/23mm, High-eyepoint	A51.3621-1023	
	WF15x/15mm, High-eyepoint	A51.3621-1515	
	WF20x/10mm, High-eyepoint	A51.3621-2010	
	WF25x/9mm, High-eyepoint	A51.3621-2509	
	WF30x/8mm, High-eyepoint	A51.3621-3008	
	WF10x/22mm, High-eyepoint, Diopter Adjustable	A51.3622-1022	
	WF15x/16mm, High-eyepoint, Diopter Adjustable	A51.3622-1516	
	WF20x/12mm, High-eyepoint, Diopter Adjustable	A51.3622-2012	
	WF25x/9mm, High-eyepoint, Diopter Adjustable	A51.3622-2509	
	WF30x/8mm, High-eyepoint, Diopter Adjustable	A51.3622-3008	
Auxiliary Lens	0.5X,W.D. 165mm	A52.3621-0.5x	
	1.5X,W.D. 45mm	A52.3621-1.5x	
	2.0X,W.D. 30mm	A52.3621-2.0x	
CCD Adapter	1.0x, For CTV	A55.3621-1.0x	
	0.50x, For 1/2" CCD	A55.3621-0.5x	
	0.35x, For 1/3" CCD	A55.3621-0.35x	
	0.50x, For 1/2" CCD	A55.3621-N0.5x	



○ Trinocular Head + Zoom Body 0.7~4.5x,
With WF10x/20mm Eyepiece



A23.3645N-B1T



○ Eyepiece Tube Adapter Included



Technical Parameters:

A23.3645N Zoom Stereo Microscope Optical Data Sheet									
Objective	Standard 1.0x		Auxiliary 0.5X		Auxiliary 1.5X		Auxiliary 2X		
W.D.	100mm		165mm		45mm		30mm		
Eyepiece	Mag	View Field	Mag	View Field	Mag	View Field	Mag	View Field	
10X/20mm	7.0X	28.6mm	3.5X	57.1mm	10.5X	19mm	14X	14.3mm	
	45.0X	4.4mm	22.5X	8.9mm	67.5X	3mm	90X	2.2mm	
15X/15mm	10.5X	21.1mm	5.25X	42.8mm	15.75X	14.3mm	21X	10.7mm	
	67.5X	3.3mm	33.75X	6.7mm	101.25X	2.2mm	135X	1.7mm	
20X/10mm	14.0X	14.3mm	7.0X	28.6mm	21.0X	9.5mm	28X	7.1mm	
	90.0X	2.2mm	45.0X	4.4mm	135.0X	1.5mm	180X	1.1mm	



Optional A52.3621 0.5x, 1.5x, 2.0x Objective

A55.3621



1.0x, For CTV



0.50x, For 1/2" CCD



0.35x, For 1/3" CCD



N0.50x, For 1/2" CCD

Optional A55.3621 CCD Adapter



Optional A55.3621 CCD Adapter

[LEARN MORE](#)



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China