



Pole Stand Stereo Zoom Microscopes Trinocular Magnification 6.7x - 45x

Our Product Introduction

for more products please visit us on cnoec.com

Basic Information

- Place of Origin: China
- Brand Name: OPTO-EDU
- Certification: CE, Rohs
- Model Number: A23.3667-B1T
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Eyepiece: WF10x/20mm
- Zoom Lens: 0.67~4.5x
- Zoom Ratio: 1:6.5
- Magnification: 6.7x-45x
- Working Distance: 105mm
- Stand: Pole Stand, No Light Source
- Highlight: **stereo zoom microscopes, stereo dissecting microscope**



More Images



Product Description

0.67-4.5x, Trinocular, WF10x/22, B1 Pole Stand + A1 Holder, Stand Size 200*255*22mm , Pole Height 240mm, Dia.32mm, No Light Source, Incl. Eyepiece Adapter



| A23.3667-B1 Zoom Stereo Microscope Specification | | B | T |
|--|---|------------|---|
| Head | Binocular Head, 45° Inclined, 360°Rotatable, ±5 Diopter Adjustable, Interpupillary Adjustment 54-75mm | O | |
| | Trinocular Head, 45° Inclined, 360°Rotatable, ±5 Diopter Adjustable, Interpupillary Adjustment 54-75mm, Eyepiece Tube Adapter | | O |
| Eyepiece | WF10x/20mm, Both with Eye-cup | | |
| Zoom Range | 0.67~4.5x | | |
| Zoom Ratio | 1:6.5 | | |
| Magnification | 6.7~45x | | |
| Working Distance | 105mm | | |
| Stand | Pole Stand, No Light Source | | |
| Optional Accessories | | Cata.No. | |
| With WF10X eyepiece; 6.7-45x continuously adjustable | | A53.3667-B | |

| | | |
|----------------|---|----------------|
| Head | Trinocular With WF10X eyepiece; 6.7-45x continuously adjustable | A53.3667-T |
| Eyepiece | WF10x/20mm, High-eyepoint | A51.3621-1020 |
| | WF10x/22mm, High-eyepoint | A51.3621-1022 |
| | WF10x/23mm, High-eyepoint | A51.3621-1023 |
| | WF15x/15mm, High-eyepoint | A51.3621-1515 |
| | WF20x/10mm, High-eyepoint | A51.3621-2010 |
| | WF25x/9mm, High-eyepoint | A51.3621-2509 |
| | WF30x/8mm, High-eyepoint | A51.3621-3008 |
| | WF10x/22mm, High-eyepoint, Diopter Adjustable | A51.3622-1022 |
| | WF15x/15mm, High-eyepoint, Diopter Adjustable | A51.3622-1515 |
| | WF20x/10mm, High-eyepoint, Diopter Adjustable | A51.3622-2010 |
| | WF25x/9mm, High-eyepoint, Diopter Adjustable | A51.3622-2509 |
| | WF30x/8mm, High-eyepoint, Diopter Adjustable | A51.3622-3008 |
| Auxiliary Lens | 0.3X, W.D. 287mm | A52.3622-0.3x |
| | 0.5X, W.D. 177mm | A52.3622-0.5x |
| | 0.7X, W.D. 113mm | A52.3622-0.7x |
| | 0.75X, W.D. 116mm | A52.3622-0.75x |
| | 1.5X, W.D. 47mm | A52.3622-1.5x |
| | 2.0X, W.D. 26mm | A52.3622-2.0x |
| CCD Adapter | 1.0x, For CTV | A55.3622-1.0x |
| | 0.50x, For 1/2" CCD | A55.3622-0.5x |



Optional A52.3622 0.3x, 0.5x, 0.7x, 0.75x, 1.5x, 2.0x Objective

A55.3621



0.50x, For 1/2" CCD



1.0x, For CTV

Optional A55.3622 CCD Adapter

A23.3667 Zoom Stereo Microscope Optical Data Sheet

| Objective | Standard 1.0x | | Auxiliary 0.5X | | Auxiliary 1.5X | | Auxiliary 2X | |
|-----------|---------------|------------|----------------|------------|----------------|------------|--------------|------------|
| W.D. | 100mm | | 165mm | | 45mm | | 30mm | |
| Eyepiece | Mag | View Field | Mag | View Field | Mag | View Field | Mag | View Field |
| 10X/20mm | 7.0X | 28.6mm | 3.5X | 57.1mm | 10.5X | 19mm | 14X | 14.3mm |
| | 45.0X | 4.4mm | 22.5X | 8.9mm | 67.5X | 3mm | 90X | 2.2mm |
| 15X/15mm | 10.5X | 21.1mm | 5.25X | 42.8mm | 15.75X | 14.3mm | 21X | 10.7mm |
| | 67.5X | 3.3mm | 33.75X | 6.7mm | 101.25X | 2.2mm | 135X | 1.7mm |
| 20X/10mm | 14.0X | 14.3mm | 7.0X | 28.6mm | 21.0X | 9.5mm | 28X | 7.1mm |
| | 90.0X | 2.2mm | 45.0X | 4.4mm | 135.0X | 1.5mm | 180X | 1.1mm |



Opto-Edu (Beijing) Co., Ltd.



0086 13911110627



sale@optoedu.com



cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China