



105mm Working Distance Stereo Optical Microscope With 0.67 - 4.5x Zoom Lens

Our Product Introduction

for more products please visit us on cnoec.com

Basic Information

- Place of Origin: China
- Brand Name: OPTO-EDU
- Certification: CE, Rohs
- Model Number: A23.3667-B1
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Eyepiece: WF10x/20mm
- Zoom Lens: 0.67~4.5x
- Zoom Ratio: 1:6.5
- Magnification: 6.7x-45x
- Working Distance: 105mm
- Stand: Pole Stand, No Light Source
- Highlight: **stereo dissecting microscope, stereo zoom microscopes**



More Images



Product Description

A23.3667-B1 Zoom Stereo Microscope Specification		B	T
Head	Binocular Head, 45° Inclined, 360°Rotatable, ±5 Diopter Adjustable, Interpupillary Adjustment 54-75mm	O	
	Trinocular Head, 45° Inclined, 360°Rotatable, ±5 Diopter Adjustable, Interpupillary Adjustment 54-75mm, Eyepiece Tube Adapter		O
Eyepiece	WF10x/20mm,Both with Eye-cup		
Zoom Range	0.67~4.5x		
Zoom Ratio	1:6.5		
Magnification	6.7~45x		
Working Distance	105mm		
Stand	Pole Stand, No Light Source		
Optional Accessories		Cata.No.	
Head	With WF10X eyepiece; 6.7-45x continuously adjustable	A53.3667-B	
	Trinocular With WF10X eyepiece; 6.7-45x continuously adjustable	A53.3667-T	
Eyepiece	WF10x/20mm, High-eyepoint	A51.3621-1020	
	WF10x/22mm, High-eyepoint	A51.3621-1022	
	WF10x/23mm, High-eyepoint	A51.3621-1023	
	WF15x/15mm, High-eyepoint	A51.3621-1515	
	WF20x/10mm, High-eyepoint	A51.3621-2010	
	WF25x/9mm, High-eyepoint	A51.3621-2509	
	WF30x/8mm, High-eyepoint	A51.3621-3008	
	WF10x/22mm, High-eyepoint, Diopter Adjustable	A51.3622-1022	
	WF15x/15mm, High-eyepoint, Diopter Adjustable	A51.3622-1515	
	WF20x/10mm, High-eyepoint, Diopter Adjustable	A51.3622-2010	
	WF25x/9mm, High-eyepoint, Diopter Adjustable	A51.3622-2509	
	WF30x/8mm, High-eyepoint, Diopter Adjustable	A51.3622-3008	
Auxiliary Lens	0.3X,W.D. 287mm	A52.3622-0.3x	
	0.5X,W.D. 177mm	A52.3622-0.5x	
	0.7X,W.D. 113mm	A52.3622-0.7x	
	0.75X,W.D. 116mm	A52.3622-0.75x	
	1.5X,W.D. 47mm	A52.3622-1.5x	
	2.0X,W.D. 26mm	A52.3622-2.0x	
CCD Adapter	1.0x, For CTV	A55.3622-1.0x	
	0.50x, For 1/2" CCD	A55.3622-0.5x	


OPTO-EDU

Zoom Stereo Microscope

A23.3667-B1

 0.67x-4.5x Zoom Body

 WF10x/20mm Eyepiece High Eyepoint

 Pole Stand Without Light Source



A23.3667 Zoom Stereo Microscope Optical Data Sheet

Objective	Standard 1.0x		Auxiliary 0.5X		Auxiliary 1.5X		Auxiliary 2X	
W.D.	100mm		165mm		45mm		30mm	
Eyepiece	Mag	View Field	Mag	View Field	Mag	View Field	Mag	View Field
10X/20mm	7.0X	28.6mm	3.5X	57.1mm	10.5X	19mm	14X	14.3mm
	45.0X	4.4mm	22.5X	8.9mm	67.5X	3mm	90X	2.2mm
15X/15mm	10.5X	21.1mm	5.25X	42.8mm	15.75X	14.3mm	21X	10.7mm
	67.5X	3.3mm	33.75X	6.7mm	101.25X	2.2mm	135X	1.7mm
20X/10mm	14.0X	14.3mm	7.0X	28.6mm	21.0X	9.5mm	28X	7.1mm
	90.0X	2.2mm	45.0X	4.4mm	135.0X	1.5mm	180X	1.1mm

PRODUCT DISPLAY



○ Eyepiece WF10X/20mm

○ 0.67x-4.5x, Magnification Up To 3.35x-270x With Optional Eyepieces And Auxiliary Objectives



A54.3640-A1






A54.3630-B1



○ B1 Pole Stand + A1 Holder



 0086 13911110627  sale@optoedu.com  cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China