



## 10x High Eyepiont Portable Stereo Microscope With 100mm Working Distance

### Our Product Introduction

#### Basic Information

- Place of Origin: China
- Brand Name: OPTO-EDU
- Certification: CE, Rohs
- Model Number: A22.3660-B1
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



#### Product Specification

- Eyepiece: WF10x/20mm
- Head: 45°Binocular Head
- Objective: Turret Objective 1x3x Or 2X/4X
- Magnification Range: 10X-80X
- Working Distance: 100mm
- Stand: Pole Stand, No Light Source
- Highlight: **stereo dissecting microscope, stereo zoom microscopes**



## Product Description

Stereo microscope is widely used in materials such as macroscopic surface observation, failure analysis, fracture analysis and other industries field

| A22.3660-B1 Stereo Microscope Specification |                                              |
|---------------------------------------------|----------------------------------------------|
| Head                                        | 45° Binocular head                           |
| Eyepiece                                    | Wide Field, High-eyepoint Eyepiece-WF10X20mm |
| Objective                                   | Turret Objective 1x3x or 2X/4X               |
| Magnification Range                         | 10X-80X                                      |
| Inter-Pupillary Distance                    | 55-75mm                                      |
| Diopter Range                               | 5mm                                          |
| Working Distance                            | 100mm                                        |
| Stand                                       | Pole Stand, No Light Source                  |



| Objective Optics Parameter |               |                 |                 |                 |                  |
|----------------------------|---------------|-----------------|-----------------|-----------------|------------------|
| Model                      | Magnification | Eyepiece        |                 |                 | Working Distance |
|                            |               | SZMEWh<br>10X20 | SZMEWh<br>15X15 | SZMEWh<br>20X10 |                  |
| A22.3660-B1-12             | 1X            | 10X20mm         | 15X15mm         | 20X10mm         | 100mm            |
|                            | 2X            | 10X20mm         | 30X7.5mm        | 40X5mm          |                  |
| A22.3660-B1-13             | 1X            | 10X20mm         | 15X15mm         | 20X10mm         |                  |
|                            | 3X            | 30X6.7mm        | 45X5mm          | 60X3.3mm        |                  |
| A22.3660-B1-24             | 2X            | 10X20mm         | 30X7.5mm        | 40X5mm          |                  |
|                            | 4X            | 40X5mm          | 60X3.75mm       | 80X2.5mm        |                  |



○ Binocular Head + Body Magnification Up  
To 10x.80x, With WF10x/20mm Eyepiece



**Opto-Edu (Beijing) Co., Ltd.**

☎ 0086 13911110627

✉ [sale@optoedu.com](mailto:sale@optoedu.com)

🌐 [cnoec.com](http://cnoec.com)

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China