



A11.0101 Laboratory Optical Microscope , Compound Light Microscope Dual Viewing Head

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: OPTO-EDU
- Model Number: A11.0101
- Minimum Order Quantity: 1
- Price: Negotiation
- Packaging Details: 1 Package/10 Sets
- Payment Terms: T/T



Product Specification

- Name: Laboratory Optical Microscope
- Type: A11.0101
- Highlight: **binocular biological microscope,
biological compound microscope**

Product Description

1

- Monocular and Dual Viewing Head Available.
- Precise finished optical system and clear image
- Steady focusing system with safty stopper
- Large smooth operating mechanical stage
- Built-in Light Source 110V or 220V

A11.0101 Specification		ABC
Magnification	40X-640X	●● -
	40X-1000X	- - ●
Head	Monocular Head, 30°inclined, 360° Rotatable	●●●
Eyepiece	Ramsden WF10X/18mm	●●●
	Achromatic 4X/0.10	●●●
Objective	Achromatic 10X/0.25	●●●
	Achromatic 40X/0.65, Spring	●●●
	Achromatic 100X/1.25, Spring, Oil	- - ●
Nosepiece	Triple Revolving Nosepiece	●● -
	Quadruple Revolving Nosepiece	- - ●
Stage	Square Stage with Specimen Clips	●● -
	Square Stage with Verinal Scale	- - ●
Focusing	Separately Coarse and Fine Focusing Knobs, with Upper Stopper	●●●
Condenser	Single Condenser, with Disc Diaphragm	● - -
	Single Condenser, with Iris Diaphragm	- ●●
Light	Mirror Dia. 48mm with Holder	●●●
Source	Base with Incandescent Lamp 220V/110V 20W	●●●
Optional Accessories		Item No.
Head	Dual Viewing Head, 360° Rotatable	TM05a
	WF16X/10mm	EP16a
Eyepiece	Ramsden WF10X/18mm with Pointer	ERq10



Opto-Edu (Beijing) Co., Ltd.



0086 13911110627



sale@optoedu.com



cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China