



## OPTO-EDU A11.0110 Biological Microscope , Monocular Compound Microscope

### Our Product Introduction

for more products please visit us on [cnoec.com](http://cnoec.com)

#### Basic Information

- Place of Origin: China
- Brand Name: OPTO-EDU
- Model Number: A11.0110
- Minimum Order Quantity: 1 pcs
- Price: USD 362.5-185.5



#### Product Specification

- Head: Binocular Or Trinocular
- Eyepiece: WF10/18mm, High Eye Point
- Objective: 4x.10x.40x.100x.
- Nosepiece: Quadruple Revolving Nosepiece
- Stage: Double Layer Mechanical Stage 130x140mm, Moving Range 30x75mm
- Condenser: N.A.1.25, Spiral Up/Down, With Iris Diaphragm
- Focusing: Coarse&Fine Focusing
- Illumination: LED Lamp
- Power: AC100-240V, 50/60Hz 0.5A
- Highlight: **binocular compound microscope, wide field microscope**



#### More Images



## Product Description

## Our Product Introduction



**Opto-Edu (Beijing) Co., Ltd.**



0086 13911110627



sale@optoedu.com



cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China

for more products please visit us on [cnoec.com](http://cnoec.com)