



40X-1600X Binocular Student Biological Microscope A11.1522-D Xsz-107bn

Our Product Introduction

for more products please visit us on cnoec.com

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A11.1522
- Minimum Order Quantity: 1 pc
- Price: Negotiation
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Total Magnification: 40X 1600X
- Eyepiece: WF10X, WF16X
- Head: Compensation Free Binocular, 45° Inclined, 360° Rotatable (White Head)
- Objectives: 4x, 10x, 40x, 100x
- Nosepiece: Quadruple
- Focusing: Coaxial Coarse And Fine Focus Adjustment, Range: 30mm
- Stage: Double Layer Mechanical Stage 142mmx132mm, Moving Range : 75mmx45mm
- Condenser: ABBE NA1.25 Condenser With Diaphragm And Filter
- Illuminator: Halogen Lamp 6V/20W Adjustment
- Highlight: **binocular biological microscope, biological compound microscope**



More Images



Product Description

Cheapest xsz-107bn 40X-1600X BInocular OPTO-EDU A11.1522-D Student Biological Microscope

Cheapest Classic Model XSZ-107BN
WF10x Eyepiece, Achromatic 4x 10x 40x100x Objectives
Double Layer Mechanical Stage 142x132mm
Halogen Lamp 6V/20W Adjustment

A11.1522 Compound Biological Microscope Model, Can Be Widely Used In Medical And Sanitary Establishments, Laboratories, Institutes, Network Of Agricultural Research, Colleges And Universities For Clinical Microscope Inspections And Used As Instructional Equipment. They Also Can Be Used In Research Work Such As Biology, Bacteriology, Cytology, Histology And Pharmaceutical Chemistry.

A 11.1522 Specification		A	B	C	D
Total Magnification	40X 1600X	•	•	•	•
Eyepiece	WF10X	•	•	•	•
	WF16X	•	•	•	•
Head	Monocular,30° Inclined,360° Rotatable	•	-	-	-
	Sliding Binocular,45° Inclined,360° Rotatable(Black Head)	-	•	-	-
	Sliding Trinocular,45° Inclined,360° Rotatable	-	-	•	-
	Compensation Free Binocular,45° Inclined,360° Rotatable(White Head)	-	-	-	•
Objectives	4x	•	•	•	•
	10x	•	•	•	•
	40x	•	•	•	•
	100x	•	•	•	•
Nosepiece	Quadruple	•	•	•	•
Focusing	Coaxial Coarse And Fine Focus Adjustment, Range:30mm	•	•	•	•
Stage	Double Layer Mechanical Stage 142mmx132mm, Moving Range : 75mmx45mm	•	•	•	•
Condenser	ABBE NA1.25 Condenser With Diaphragm And Filter	•	•	•	•
Illuminator	Halogen Lamp 6V/20W Adjustment	•	•	•	•

OPTO-EDU


Compound Biological Laboratory Microscope

A11.1522

- Cheapest Classic Model XSZ-107BN
- WF10x Eyepiece, Achromatic 4x 10x 40x100x Objectives
- Double Layer Mechanical Stage 142x132mm
- Halogen Lamp 6V/20W Adjustment

Don't Miss Best Microscope For You!

A11.1522 COMPOUND BIOLOGICAL LABORATORY MICROSCOPE



A11.1522 Compound Biological Laboratory Microscope Practical Application



A11.1522-A



A11.1522-B




A11.1522-C



A11.1522-D

A11.1522 Biological Microscope



PRODUCT DISPLAY



01 HEAD

Compensation Free Binocular, 45° inclined, 360° Rotatable (White Head)

02 OBJECTIVE

Quadruple Nosepiece, Achromatic 4x, 10x, 40x, 100x Objectives Included



03 WORKING STAGE

Double Layer Mechanical Stage:
132mmx142mm,
Moving Range : 75mmx45mm

04 ILLUMINATOR

Halogen Lamp 6V/20W Adjustment, Abbe Condenser With Iris Diaphragm, With Filter Holder, And Moving Up/Down By Adjust Knob. Center Adjustable By 2 Screws On Condenser Holder



05 FOCUSING

Coaxial Coarse And Fine Focus Adjustment,
Range: 30mm

PRODUCT DESCRIPTION





PACKING AND SHIPPING

Packing : Foam + Wooden Box + Carton, 34.5*30.5*49CM, G.W.: 10Kgs (Wooden Box Optional)

Shipping Time : 3-10 Days From Stock



Opto-Edu (Beijing) Co., Ltd.



0086 13911110627



sale@optoedu.com



cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China