



Optical 20x 40x Binocular Stereo Microscope Sunshine Student Boom Stand Inspection Electron Pcb Repair

Our Product Introduction

for more products please visit us on cnoec.com

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE,
- Model Number: A22.1237
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Head: Rotatable 360°, Inclination 45° Binocular, With Locked-in Eyepiece, Left & Right Diopter Adjustable $\pm 5^\circ$, Head Interpupillary 55-75mm
- Eyepiece: Extra Wide Field, Plan WF10x/22mm, High Eyepoint, Eyepiece Tube $\Phi 30\text{mm}$ WF10X
- Objective: Three Kind Of Objective 90° Rotatable For Selection: 1X2X, 1X3X, 2X4X
- Working Stage: $\Phi 60$ White/Black And $\Phi 60$ Frosted Glass
- Power Supply: AC 100V-240V Wide Range, CE Approval
- Light Source: No Light Source Or Reflect 1WLED, Transmit 7x20mA
- Highlight: **stereo zoom microscopes,**
stereo binocular microscopes,
stereo dissecting microscope



More Images



Product Description

OPTO-EDU Cheap Student Boom Stand Inspection Electron Pcb Repair Optical 20x 40x Binocular Stereo Microscope Sunshine

Stereo microscope is widely used in teaching demonstration, bio-engineering, IT industry testing and other fields. A22.1237 has a high resolution and large field of view clarity, so it is also good at electronics industry and precision machinery industry parts assembly and testing, agricultural seed inspection and so on.

Model	P	PA	PB	TA	TB
Head	Rotatable 360°, Inclination 45° Binocular, With Locked-in Eyepiece, Left & Right Diopter Adjustable ±5°, Head Interpupillary 55-75mm				
Eyepiece	WF10X				
Objective	Three Kind Of Objective 90° Rotatable For Selection: 1X2X, 1X3X, 2X4X				
Light Source	No Light Source	Reflect 1WLED, Transmit 7x20mA			
Stand	Pole Stand	Pole Stand Touch Switch	Pole Stand Dial Switch	Track Stand Touch Switch	Track Stand Dial Switch
Working Stage	Φ60 White/Black And Φ60 Frosted Glass				
Power Supply	AC 100V-240V Wide Range, CE Approval				
Optional Accessories					
Eyepiece	WF5X				
	SWF10X				
	WF15X				
	WF20X				

Objective Magnification	Magnification	Field Of View(mm)	Magnification	Field Of View(mm)	Magnification	Field Of View(mm)	Magnification	Field Of View(mm)	Working Distance(mm)
1X	5X	20	10X	20	15X	15	20X	10	76
2X	10X	10	20X	10	30X	7.5	40X	6.5	76
3X	15X	7	30X	6.7	45X	5	60X	4.3	76
4X	20X	5	40X	5	60X	3.7	80X	3.2	76



OPTO-EDU


Stereo Microscope

A22.1237

- Head 360° Rotatable With Locked-in Eyepiece
- WF10x/20mm High Eyepoint Eyepiece, Diopter Adjustable
- Objective 90° Rotatable: 1X2X, 1X3X, 2X4X
- Sharp Stereo Erect Image Over A Wide Field Of View

Don't Miss Best Microscope For You!

A22.1237 STEREO MICROSCOPE



 A22.1237-TA Track Stand+Touch Switch Binocular Stereo Microscope

A22.1237-P Pole Stand



A22.1237-PA Pole Stand+Touch Switch



A22.1237-PB Pole Stand+Dial Switch



A22.1237-TB Track Stand+Dial Switch



The coin under the stereo microscope:



PRODUCT DISPLAY




01 BIBOCULAR HEAD

Rotatable 360°, Inclination 45° Binocular, With Locked-in Eyepiece, Left & Right Diopter Adjustable $\pm 5^\circ$, Head Interpupillary 55-75mm

02 OBJECTIVE

Three Kind Of Objective 90° Rotatable For Selection: 1X2X, 1X3X, 2X4X



03 FOCUSING

Coarse Focus, Focusing Range: 40mm

04 LIGHT SOURCE

Reflect 1W LED Light Source



05 STAGE PLATE

Φ60 White/Black And Frosted Glass

06 TOUCH SWITCH

Light Source Brightness Adjustment



07 POWER & SWITCH KNOB

AC 100V-240V Wide Range, CE Approval



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China