



## Optical Pcb Mobile Repair Zoom Stereo Lcd Microscope 3.0M CMOS Usb2.0 Led Lighted Electronic HD Screen Display

Our Product Introduction

for more products please visit us on cnoec.com

### Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A36.1210
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



### Product Specification

- LCD Screen: 9" Full HD LCD Pixel Size 2.75μm×2.75μm Max Fps 30fps, USB2.0 SD Card Power Switch, Can Rotate 0° 100°
- Camera: 1/3" 3.0M CMOS Sensor
- Objective: Zoom Lens 1× 4× + Reduction Les 0.3×
- Working Distance: W.D. 90.5mm
- Magnification: Digital Mag. 38.1× Total Mag. 11.5× 46×. Turn Middel Wheel Of Mouse Can Digital Zoom 1× 10×
- Light Source: Upper Light 21×20mA LED Bottom Light 7×20mA LED, With Brightness Adjust Panel
- Power Supply: 100V 240V / 12V 1A or 2A USB2.0 Mouse
- Software: Mouse Connect To LCD USB2.0 Move Mouse To Left Edge Of Screen Can See Menu
- Highlight: digital trinocular microscope, usb digital microscope, infinity plan microscope



### More Images



### Product Description

#### Opto-edu Led Lighted Mobile Usb2.0 3.0m Cmos 30fps New Best Price Zoom Electronic 1080p Lcd Digital Display Video Microscope

#### Product Description

A36.1210 LCD Stereo Microscope, 9"	
LCD Screen	9" Full HD LCD, Pixel Size 2.75μm×2.75μm, Max Fps 30fps, USB2.0, SD Card, Power Switch, Can Rotate 0°~100°
Camera	1/3" 3.0M CMOS Sensor
Objective	Zoom Lens 1×~4× + Reduction Les 0.3×
Working Distance	W.D. 90.5mm
Magnification	Digital Mag. 38.1×, Total Mag. 11.5×~46×. Turn Middel Wheel Of Mouse Can Digital Zoom 1×~10×
Light Source	Upper Light 21×20mA LED, Bottom Light 7×20mA LED, With Brightness Adjust Panel
Power Supply	100V~240V / 12V 1A (or 2A), USB2.0 Mouse
Software	Mouse Connect To LCD USB2.0, Move Mouse To Left Edge Of Screen Can See Menu

Our Product Introduction



**OPTO-EDU**

## LCD Stereo Microscope, 9"

**A36.1210**

- 9" HD LCD View Angle Adjustable 0°~100°
- 3.0M CMOS Digital Camera
- Zoom Lens 1x~4x, Total Magnification 11.5x~460x
- Up/Bottom LED Light Source Brightness Adjustable
- USB Mouse Control + Measuring Function + SD Card Slot



**Don't Miss Best Microscope For You!**

### A36.1210 LCD STEREO MICROSCOPE



### PRODUCT DISPLAY



#### 01 LCD SCREEN

9" Full HD LCD, Pixel Size 2.75μm×2.75μm, Max Fps 30fps, 1/3" 3.0M CMOS Sensor, USB2.0, SD Card, Power Switch, Can Rotate 0°~100°

#### 02 OBJECTIVE

Zoom Lens 1x~4x + Reduction Les 0.3x







### 03 FOCUSING

Coarse Focus, Working Distance: 90.5mm

### 04 LIGHT SOURCE

Upper Light 21×20mA LED, With Brightness Adjust Panel



### 05 WORKING STAGE

Track Stand Working Stage, Bottom Light 7×20mA LED, With Brightness Adjust Panel

### 06 SWITCH

Power Switch And Touch Switch To Adjust Brightness Of Both Up/Bottom Light Source



### 07 POWER SUPPLY

100V~240V / 12V 1A (or 2A) 、USB2.0 Mouse



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China