



## Handheld Mini Digital Metallurgical Microscope OPTO-EDU A13.2501-B 100 - 400x

Our Product Introduction

for more products please visit us on cnoec.com

### Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE,
- Model Number: A13.2501
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



### Product Specification

- Total Magnification: 100-400x
- Head: Monocular
- Eyepiece: WF10X
- LWD Plan Objective (no Cover Glass): PLL 10X/0.25 WD:7.3mm & PLL 40X/0.65 WD:0.5mm
- Mechanical Tube Length: 160mm
- Focusing: Coaxial Coarse/fine Focus System, Range:36mm
- Light Source: Coaxial Illumination, Build-in LED Rechargeable, Brightness Adjustable
- Packing: Portable Aluminium Case
- Highlight: reflected light microscope, halogen lamp microscope



### More Images



### Product Description

#### OPTO-EDU A13.2501-B 100-400x Handheld Mini Metallurgical Portable Metallurgical Microscope

##### Product Description

Small & Simple, with bright white LED light built-in, and high quality image!

A13.2501-B Portable Metallurgical Microscope		
Total Magnification	100-400x	
Head	Monocular	
Eyepiece	WF10X	
LWD Plan Objective	PLL 10X/0.25 WD:7.3mm	
(no cover glass)	PLL 40X/0.65 WD:0.5mm	
Mechanical Tube Length	160mm	
Focusing	Coaxial coarse/fine focus system, range:36mm	
Light Source	Coaxial illumination, build-in LED rechargeable, brightness adjustable	
Packing	Portable aluminium case	
A13.2501 Portable Metallurgical Microscope Optional Accessories		
Adapter	For Canon SLR Digital Camera	A55.2502
Stand	Magnet Stand	A54.2501
	PLL 20x	A5M.0234-20
	PLL 60x	A5M.0234-60
	PLL 80x	A5M.0234-80
	PLL 100x	A5M.0234-100
LWD Plan Objective		
(no cover glass)		

Our Product Introduction



LWD Plan Objective



Focusing



Portable Aluminum Case



A55.2502 Canon SLR Camera Adapter





Optional Part: A54.2501 Magnet Stand





Company Information





OPTO-EDU is you're looking for, isn't it?

OPTO-EDU CNOPTec



Our Team And Our Factory.

OPTO-EDU, as one of the most professional manufacturer and supplier of microscope in China, our sub-brand CNOPTec series high end biological, laboratory, polarizing, metallurgical, fluorescence microscopes, CNCOMPARISON series forensic microscope, A63 series SEM microscope, and .49 series digital camera, LCD camera are very popular in the world market.



I'm OPTO-EDU, nice to meet you!



Trust us, we will give you a guarantee.

Direct factory price!  
3 years quality warranty!  
Payment by T/T, PayPal, West Union to LC!  
Ranked No.1 supplier for microscope on alibaba.com!  
Special price support OPTO-EDU/CNOPTec distributor!

Opto-Edu (Beijing) Co.,Ltd.

Transaction Level

41 Transaction(s)(6months)

\$ 1,400,000+

Response Rate

100.0%

The Supplier Supports Trade Assurance-A Free Service That Protects Your Orders From Payment To Delivery

Supplier's Trade Assurance Limit: US \$379,000






**Opto-Edu (Beijing) Co., Ltd.**

 0086 13911110627

 [sale@optoedu.com](mailto:sale@optoedu.com)

 [cnoec.com](http://cnoec.com)

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China