



A43.1901 C Mount 1/3" Ccd Surgical Operating Microscope 2.9x - 21.7x Colposcope

Our Product Introduction

for more products please visit us on cnoec.com

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE,
- Model Number: A43.1901
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Total Magnification: 2.9x-21.7x
- Magnification Changer: 5-step Magnification Change
- Beamsplitter: 50: 50 Beamsplitter
- TV Adapter: C-mount 1/3" CCD Attachable
- Highlight: **colposcope surgical operating microscope, 21.7x surgical operating microscope, 2.9x surgical microscope**



More Images



Product Description

Specifications					
Model	A43.1901-A	A43.1901-B	A43.1901-C	A43.1901-D	A43.1901-E
Observation angle	45° Tilted Binocular Observation Unit				
Magnification of eyepiece	16X				
Focal length of objective	F=230mm & 280mm (Connecting screw: M58x0.75)			F=280mm (Connecting screw: M57x0.75)	
Working distance	215mm & 260mm			265mm	
Magnification changer	0.4X, 0.66X, 1X, 1.5X, 2.5X 5-step magnification change			No	
Total magnification	2.9X-21.7X			11X	
Linear field	11mm - 84mm			24mm	
Beamsplitter	50:50 Beamsplitter			No	
Demonstrator	Total magnification:2.1X-16.4X			Optional	
TV adapter	C-mount 1/3" CCD attachable				No
Photo adapter	Optional				No
Fine focusing distance	20mm	60mm			20mm
Balancing arm	Counterweight adjustable	No	Counterweight adjustable		
Illumination system	Coaxial illumination with cold light source fiber-optic cable 12V 100W halogen lamp, brightness adjustable				
Filter	Green filter with frame, to be fitted on the objective when it is used				
Fiber-optic cable	1300mm				
Stand	short column fitted on the round base with 5 jaws and a round counterweight ring	Column height adjustable fitted on the round base with 3 jaws			short column fitted on the round base with 5 jaws



A43.1901-A Colposcop

www.cnoec.com



A43.1901-B Colposcop



A43.1901-A+Demonstrator+Camera+Computer



A43.1901-A+Camera Adapter+Camera+Computer



A43.1901-B+Demonstrator



A43.1901-B+Demonstrator+Computer

A43.1901 Colposcop



Opto-Edu (Beijing) Co., Ltd.



0086 13911110627



sale@optoedu.com



cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China