



A24.1201 40x Stereo Jewelry Microscope / Gem Microscope Dark Field Halogen Lamp

Our Product Introduction

for more products please visit us on cnoec.com

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A24.1201
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Head: Inclined Binocular Head, 360° Rotatable
- Eyepiece: WF10x
- Objective: 2x4x
- Incident Light: Halogen Lamp 12V/10W, Brightness Adjustable
- Transmitted Light: Halogen Lamp 12V/10W, Brightness Adjustable
- Working Stage: Dark Field Working Stage, With Diaphragm & Jewel Tweezer
- Highlight: **laboratory binocular microscope, gemological microscope**



More Images




Product Description

40x Stereo Jewelry Microscope / Gem Microscope Dark Field Halogen Lamp A24.1201

Gem microscope is used in jewelry identification to observe the intrinsic and surface characteristics of the sample and is invisible to the naked eye. It is an indispensable gem identification instrument for the identification of synthetic products, optimized treatments and artifacts. It also has extremely important applications in diamond grading, such as observing clarity characteristics, observing and evaluating the color grade of diamonds.

A24.1201 Jewery Microscope	
Head	Inclined Binocular Head,360° Rotatable
Eyeiece	WF10x
Objective	2x4x
Incident Light	Halogen Lamp 12V/10W, Brightness Adjustable
Transmitted Light	Halogen Lamp 12V/10W, Brightness Adjustable
Working Stage	Jewel Tweezer Dark Field Working Stage, With Diaphram
A24.1201 Jewery Microscope , Optional Accessories	
Eyeiece	WF5x E51.1221-05
	WF15x E51.1221-15
	WF20x E51.1221-20

Company Information



📍 OPTO-EDU is you're looking for, isn't it?

OPTO-EDU *cnoptec*

<p>25+ Yeas Professional Experience Know Micro- scopes Best</p> <p>150+ 150+ Micro- scope & Accessory Manufactur- ers From</p> <p>25+</p>	<p>200+ Hot Sale & Newest Update Mic- scope Model Every Month</p> <p>720+ 720+ Cus- tomers From All Over World Keep Rising</p> <p>200+</p>	<p>1200+ Micro- scope Models To Creat One-Stop Purchase Platform</p> <p>1200+</p>	<p> Top Seller For Micro- scope on Alibaba Recommand</p>	<p>OPTO-EDU As Famous Brand For Mi- croscope In China & USA</p> <p> Most Professional Sale Team Supply Best Microscope</p>
---	--	--	--	--

Our Team And Our Factory.

OPTO-EDU, as one of the most professional manufacturer and supplier of microscope in China, our sub-brand CNOPTEC series high end biological, laboratory, polarizing, metallurgical, fluorescene microscopes, CNCOMPARISON series forensic microscope, A63 series SEM microscope, and .49 series digital camera, LCD camera are very popular in the world market.



😊 I'm OPTO-EDU, nice to meet you!



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China