



40x - 1000x Biological Microscope Opto-Edu A11.1123 WF10X With Pointer

Our Product Introduction

for more products please visit us on cnoec.com

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A11.1123
- Minimum Order Quantity: 1 pc
- Price: Negotiation
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Head: Sliding Binocular Head, 45°Inclined, 360°Rotatable
- Eyepiece: WF10X With Pointer
- Objective: Achromatic 4x, 10x, 40x, 100x
- Nosepiece: Quadruple Revolving Nosepiece
- Stage: Double Layer Mechanical Stage ,Size 100x90mm
- Condenser: NA1.25 Condensor With Iris Diaphragm
- Focusing: Coaxial Coarse And Fine Focusing Knobs
- Highlight: **high powered microscope, binocular biological microscope, biological compound microscope**



More Images



Product Description

Opto-Edu A11.1123 40x-1000x Economical Binocular Optical Microscope

Product Description

A 11.1123 Economical Binocular Biological microscope	
This scope is an ideal school microscope, it is with LED lighting and fully portable with no need for messy classroom wires! Arm with handle hole is easy to carry	
A11.1123 Economical Binocular Biological Microscope Specification	
Head	Sliding Binocular Head, 45°Inclined, 360°Rotatable
Eyeiece	WF10X with pointer
Objective	Achromatic 4X
	Achromatic 10X
	Achromatic 40X, Spring
	Achromatic 100X, Spring/Oil
Nosepiece	Quadruple Revolving Nosepiece
Stage	Double layer Mechanical Stage ,Size 100x90mm
Condenser	NA1.25 Condensor with Iris Diaphragm
Focusing	Coaxial Coarse and Fine Focusing Knobs
Light Source	Top and Bottom LED light with rechargeable function (if shipped by air, battery will not be included)
	Power supply adapter



A11.1123 Biological Microscope




A11.1123 Side View

A11.1123 Back View



Upper LED Light Source



Coaxial Coarse & Fine Focusing Adjust




Bottom LED Light Source




Both Eyepiece Tube Diopter Adjustable

Company Information



📍 OPTO-EDU is you're looking for, isn't it?

OPTO-EDU *cnoptec*



Our Team And Our Factory.

OPTO-EDU, as one of the most professional manufacturer and supplier of microscope in China, our sub-brand CNOPTec series high end biological, laboratory, polarizing, metallurgical, fluorescence microscopes, CNCOMPARISON series forensic microscope, A63 series SEM microscope, and .49 series digital camera, LCD camera are very popular in the world market.



😊 I'm OPTO-EDU, nice to meet you!



Opto-Edu (Beijing) Co., Ltd.



0086 13911110627



sale@optoedu.com



cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China