



1.3 Mega Pixels Digital Optical Microscope 20x - 400x A34.5007 1.3M USB2.0 Digital Microscope

Our Product Introduction

for more products please visit us on cnoec.com

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A34.5007
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Image Sensor: 1.3 Mega Pixels, (Interpolated To 2.0M)
- Still Capture Resolution: 1600*1200(2M), 1280*1024, 1280*960, 1024*768, 800*600, 640*480, 352*288, 320*240
- Video Capture Resolution: 1600*1200(2M), 1280*1024, 1280*960, 1024*768, 800*600, 640*480, 352*288, 320*240
- Focus Range: Manual Focus From 10mm To 500mm
- Frame Rate: Max 30fps Under 600 Lus Brightness
- Magnification Ratio: 20x~400x
- PC Interface: USB2.0
- Video Format: AVI
- Highlight: digital usb microscope, portable digital microscope



More Images



Product Description

A34.5007 20x-400x 1.3M USB digital microscope

Product Description

Digital Microscope, 400X, 1.3M

--Magnification 20X~400X


--Resolution 1.3M, Upto 2.0M By Software

--8 LED, Adjustable By Control Wheel

--Manual Focus From 10mm To 500mm

A34.5007 USB Digital Microscope	
Image Sensor	1.3 Mega Pixels, (Interpolated to 2.0M)
Still Capture Resolution	1600*1200(2M), 1280*1024, 1280*960, 1024*768, 800*600, 640*480, 352*288, 320*240
Video Capture Resolution	1600*1200(2M), 1280*1024, 1280*960, 1024*768, 800*600, 640*480, 352*288, 320*240
Focus Range	Manual Focus From 10mm To 500mm
Frame Rate	Max 30fps Under 600 Lus Brightness
Magnification Ratio	20x~400x
Video Format	AVI
Photo Format	JPEG / BMP
Light Source	8 LED, Adjustable By Control Wheel
PC Interface	USB2.0
Power Source	5V DC From USB Port
Operation System:	Windows 2000/XP/Vista/Win7 Mac OS X 10.5 or Above
OSD Language	English, German, Spanish, Korean, French, Russian
Bundle Software	MicroCapture With Measurement
Product Size	125mm (L) x 33mm (R)
Product Weight	160g
Packing List	Microscope Metal Stand Software CD User's Manual
Package	20Pcs/Carton, 10Kgs/Carton, 46(L)*35(W)*30.0(H)cm
Certificate	With CE, FCC Approval, RoHS Compliance









A34.5007 USB

Packaging & Shipping

Packaging:




Foam + Carton

Shipping:

By air, ship or express.




Company Information



📍 OPTO-EDU is you're looking for, isn't it?

OPTO-EDU *cnoptec*



Our Team And Our Factory.

OPTO-EDU, as one of the most professional manufacturer and supplier of microscope in China, our sub-brand CNOPTec series high end biological, laboratory, polarizing, metallurgical, fluorescence microscopes, CNCOMPARISON series forensic microscope, A63 series SEM microscope, and .49 series digital camera, LCD camera are very popular in the world market.



😊 I'm OPTO-EDU, nice to meet you!



Opto-Edu (Beijing) Co., Ltd.



0086 13911110627



sale@optoedu.com



cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China