

30fps A59.2302 Microscope Accessories 1.3M / 2M Mini Handheld Usb Digital Microscope

Our Product Introduction

for more products please visit us on cnoec.com

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A59.2302
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- USB Device: USB 2.0
- Magnification: 10x---30x And 200x (17" Monitor)
- Frame Rate: 30fps At 640*480
- Sensitivity: 1.0V/lux-sec@550nm
- Spectrum Response: 400nm---1000nm
- Shutter Type: ERS(electronic Rolling Shutter)
- Highlight: microscope accessory, parts of the microscope, compound microscope parts



More Images



Product Description

OPTO-EDU A59.2302 1.3M / 2M USB Mini Handheld Usb Digital Microscope

Product Description

This digital microscope is a hand-held digital magnifier. By connected with a PC through a USB2.0 port, the images and videos of the specimen are shown on the screen, and captured, edited by the software attached. With the special design of adjustable magnification, the fine target can be found easy in low magnification and its detailed information can be inspected in high magnification. The advanced design in IC project enables This scope run without a driver, just like an insert-plug movable hard-disk. The newly developed software attached supports the working model of video recording in high resolution (1.3MP or 2.0MP). This scope is an extremely useful tool for parts inspection and precision testing which is widely used in machinery, electronics, printing, commodity inspection, quality control and in areas like textiles industry, security, health, education and family life



| | | |
|-------------------|---|------------------------------|
| Sepcification | A59.2302-A | A59.2302-B |
| Sensor Parameter | 1/4" 1.3 Mega Pixels CMOS | 1/4" 2.0Mega Pixels CMOS |
| USB Device | USB 2.0 | |
| Magnification | 10x---30x and 200x (17" Monitor) | |
| Frame Rate | 30fps at 640*480 | |
| Preview Model | 1280*1024, 1024*768,640*480, 320*240,160*120 | 1600*1200,1280*1024, 640*480 |
| Max Resolution | 1280*1024 | 1600*1200 |
| Sensitivity | 1.0V/lux-sec@550nm | |
| Magnification | 10*---30* and 200* (17" Monitor) | |
| Spectrum response | 400nm---1000nm | |
| Shutter Type | ERS(electronic rolling shutter) | |
| Exposure | Auto | |
| White Balance | Auto | |
| USB cable length | 1.5m | |
| Light | Built-in 8pcs LEDs | |
| Video Format | AVI | |
| Operating System | Windows--XP2,Vista | |
| Basic PC | Pentium4 2.0GHz or more 512MB RAM or more USB 2.0 Port 17" Monitor or more big | |



USB 2.0 interface



Package




Built-in 8pcs LEDs




Package

Company Information



📍 OPTO-EDU is you're looking for, isn't it?

OPTO-EDU *cnoptec*



Our Team And Our Factory.

OPTO-EDU, as one of the most professional manufacturer and supplier of microscope in China, our sub-brand CNOPTec series high end biological, laboratory, polarizing, metallurgical, fluorescence microscopes, CNCOMPARISON series forensic microscope, A63 series SEM microscope, and .49 series digital camera, LCD camera are very popular in the world market.



😊 I'm OPTO-EDU, nice to meet you!



👍 Trust us, we will give you a guarantee.

Direct factory price!
3 years quality warranty!
Payment by T/T, PayPal, West Union to LC!
Ranked No.1 supplier for microscope on alibaba.com!
Special price support OPTO-EDU/CNOPTEC distributor!

Opto-Edu (Beijing) Co., Ltd.

Transaction Level

41 Transaction(s)(6months) Response Rate
\$ 1,400,000+ 100.0%

The Supplier Supports Trade Assurance-A Free Service That
 Protects Your Orders From Payment To Delivery.
 Supplier's Trade Assurance Limit: US \$379,000



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627 ✉ sale@optoedu.com 🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China