



## 0.2x - 4.9x 3.0M Compound Optical Microscope OPTO-EDU A18.1814 Multi Function Comparator

Our Product Introduction

for more products please visit us on [cnoec.com](http://cnoec.com)

### Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A18.1814
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



### Product Specification

- OPTO-EDU A18.1814 Multi 0.2x~4.9x Function Comparator:
- Working Distance: 75~90mm
- View Field: Dia. 1.5~60mm
- CCD: Built-in Digital CCD Built-in, 3.0M, Video/USB Output
- Software: Professional Forensic Image Analysis Software Crime Image
- PC: 14.1" TFT Notebook, With Software Installed
- Highlight: **binocular compound microscope, wide field microscope**



### More Images



## Product Description

### 0.2x~4.9x 3.0M With Software OPTO-EDU A18.1814 Multi Function Comparator

A18.1814 GSB-4 Multi Function Comparator Specification	
Magnification	0.2x~4.9x
Objective	Zoom objective 0.7x~4.9x, zoom ratio 1:7 Auxiliary lens 0.3x
Working Distance	75~90mm
View Field	Dia.1.5~60mm
Bullet Holder	Small bullet holder set Motorized bullet holder set
Working Stage	Titatable working stage, for bullet and small object Large motorized working stage, for document and large object Flourscence transmit light stage
Light Source	Excite & receive light source wave length: 254nm (Ultraviolet), 385nm (Ultraviolet), 540nm (Yellow), 680nm (Red), 800nm (Infrared), 960nm (Infrared), 1100nm (Infrared) Cold light source, for temperature sensitive objects viewing Holegen light source Flourscence ring light source Flourscence transmit light stage
CCD	Built-in digital CCD built-in, 3.0M, video/USB output
Software	Professional forensic image analysis software Crime Image
PC	14.1" TFT notebook, with software installed
A18.1814 GSB-4 Multi Function Comparator Packing List	
Main Body	1 set
Auxiliary lens 0.3x	1 pc
Magnification correcting lens	1 pc
Small bullet holder set	1 set
Motorized bullet holder set	1 set
Titatable working stage, for bullet and small object	1 set
Large motorized working stage, for document and large object	1 set
Flourscence transmit light stage	1 set
Flourscence ring light source	1 set
Video cable	1 pc
Power cord	1 pc
Holgen lamp	2 pcs
2A Fuse	2 pcs
Lens cleaning paper & Brush	1 set
Instruction Manual	1 book





**Opto-Edu (Beijing) Co., Ltd.**

☎ 0086 13911110627

✉ [sale@optoedu.com](mailto:sale@optoedu.com)

🌐 [cnoec.com](http://cnoec.com)

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China