



Digital Stereo Microscope

Our Product Introduction

for more products please visit us on cnoec.com

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE,
- Model Number: A32.3602-STL2
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Head: Binocular & Trinocular Head, 45° Inclined, 360°Rotatable, ± 5 Diopter Adjustable, Interpupillary Adjustment 54-75mm
- Eyepiece: WF10x/20mm
- Zoom Lens: 0.7~4.5x
- Zoom Ratio: '1:6.5
- Magnification: 7x-45x, Up To 3.5x-180x With Optional Eyepieces And Auxiliary Objectives
- Working Distance: 100mm, With Auxiliary Objectives 30~165mm
- Highlight: **digital usb microscope,**
portable digital microscope



More Images



Product Description

A32.3602-STL2 Metallurgical Optical Microscope / Digital Stereo Microscope

Stereo Microscope Magnification 7-45x, Can Up To 3.5-180x

Camera 3.5M, 1080p 30 fps, Dual Output USB/HDMI

Fine Objective Lens Technology, The Observation Effect Is More Real

Special Manual Focus On The Top Of Camera To Extend Working Distance


Dual Arm Boom Stand, More Flexible To Use



OPTO-EDU

A32.3602-STL2
Digital Stereo Microscope
HDMI + USB

- Stereo Microscope Magnification 7-45x, Can Up To 3.5-180x
- Camera 3.5M, 1080p 30 fps, Dual Output USB/HDMI
- Fine Objective Lens Technology, The Observation Effect Is More Real
- Special Manual Focus On The Top Of Camera To Extend Working Distance
- Dual Arm Boom Stand, More Flexible To Use



Optional A51.3621 Wide Field Eyepiece



WF10X/22mm



WF15X/15mm



WF25X/9mm



WF10X/23mm



WF20X/10mm

Optional A51.3622 Wide Field Eyepiece, Diopter Adjustable



WF10X/22mm



WF30X/8mm



WF25X/9mm



WF15X/16mm



WF20X/12mm

Optional A52.3621 0.5x, 1.5x, 2x Objective



A23.3645 Zoom Stereo Microscope Optical Data Sheet									
Objective	Standard 1.0x		Auxiliary 0.5X		Auxiliary 1.5X		Auxiliary 2X		
W.D.	100mm		165mm		45mm		30mm		
Eyepiece	Mag	View Field	Mag	View Field	Mag	View Field	Mag	View Field	
10X/20mm	7.0X	28.6mm	3.5X	57.1mm	10.5X	19mm	14X	14.3mm	
	45.0X	4.4mm	22.5X	8.9mm	67.5X	3mm	90X	2.2mm	
15X/15mm	10.5X	21.1mm	5.25X	42.8mm	15.75X	14.3mm	21X	10.7mm	
	67.5X	3.3mm	33.75X	6.7mm	101.25X	2.2mm	135X	1.7mm	
20X/10mm	14.0X	14.3mm	7.0X	28.6mm	21.0X	9.5mm	28X	7.1mm	
	90.0X	2.2mm	45.0X	4.4mm	135.0X	1.5mm	180X	1.1mm	

Optional A55.3621 CCD Adapter



A23.3645-STL2 Zoom Stereo Microscope Specification		B	T
Head	Binocular Head, 45° Inclined, 360°Rotatable, ±5 Diopter Adjustable, Interpupillary Adjustment 54-75mm	O	
	Trinocular Head, 45° Inclined, 360°Rotatable, ±5 Diopter Adjustable, Interpupillary Adjustment 54-75mm, Eyepiece Tube Adapter		O
Eyepiece	WF10x/20mm		
Zoom Lens	0.7~4.5x		
Zoom Ratio	1:6.5		
Magnification	7x-45x, Up To 3.5x-180x With Optional Eyepieces And Auxiliary Objectives		
Working Distance	100mm, With Auxiliary Objectives 30~165mm		
Field Of View	1~95.2mm		
Stand	Dual Arm Boom Stand, No Light Source		
Optional Accessories		Cata.No.	
Head	Binocular Head + WF10x/20mm	A53.3645-B	
	Trinocular Head + WF10x/20mm	A53.3645-T	
Eyepiece	WF10x/20mm, High-eyepoint	A51.3621-1020	
	WF10x/22mm, High-eyepoint	A51.3621-1022	
	WF10x/23mm, High-eyepoint	A51.3621-1023	
	WF15x/15mm, High-eyepoint	A51.3621-1515	
	WF20x/10mm, High-eyepoint	A51.3621-2010	
	WF25x/9mm, High-eyepoint	A51.3621-2509	
	WF30x/8mm, High-eyepoint	A51.3621-3008	
	WF10x/22mm, High-eyepoint, Diopter Adjustable	A51.3622-1022	
	WF15x/16mm, High-eyepoint, Diopter Adjustable	A51.3622-1515	
	WF20x/12mm, High-eyepoint, Diopter Adjustable	A51.3622-2010	
	WF25x/9mm, High-eyepoint, Diopter Adjustable	A51.3622-2509	
	WF30x/8mm, High-eyepoint, Diopter Adjustable	A51.3622-3008	
Auxiliary Lens	0.5X, W.D. 165mm	A52.3621-0.5x	
	1.5X, W.D. 45mm	A52.3621-1.5x	
	2.0X, W.D. 30mm	A52.3621-2.0x	
CCD Adapter	1.0x, For CTV	A55.3621-1.0x	
	0.50x, For 1/2" CCD	A55.3621-0.5x	
	0.35x, For 1/3" CCD	A55.3621-0.35x	
	0.50x, For 1/2" CCD	A55.3621-N0.5x	



Opto-Edu (Beijing) Co., Ltd.



0086 13911110627



sale@optoedu.com



cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China

