



XSZ-107BN Student Biological Microscope

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A11.107
- Minimum Order Quantity: 1 pc
- Price: Negotiation
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Applications: School|Laboratory
- Certification: CE|Rohs
- Drawtube: Binocular|Monocular|Trinocular
- Max Magnification: 1600x
- Product Name: Microscope
- Theory: Biological Microscope
- Highlight: **high powered microscope ,
binocular biological microscope,
biological compound microscope**



Product Description

Wide Field Biological Microscope 1600x With Plan - Concave Mirror A11.1007

Description:

- 195mm Achromatic Objectives, Mechanical Tube Length 160mm
- WF10x/18mm, WF16x/11mm Wide Field Eyepiece
- Monocular, Binocular, Trinocular Head Optional
- Low Position Coaxial Coarse & Fine Focus System
- Double Layer Mechanical Stage

Specifications:

A11.1007 Biological Microscope		16	27	17	17E	17C
Magnification	40x~1600x					
Head	Monocular Head, 45° Inclined, 360° Rotatable	o				
	Seidentopf Binocular Head, 45° Inclined, 360° Rotatable		o			
	Sliding Binocular Head, 45° Inclined, 360° Rotatable			o		
	Sliding Binocular Head, 45° Inclined, With Photo Adapter				o	
	Sliding Binocular Head, 45° Inclined, With C-Mount					o
Eyepiece	Wide Field WF10x/18mm					
	Wide Field WF16x/11mm					
Nosepiece	Quadruple Nosepiece					
Objective	Achromatic 4x/0.10, 10x/0.25, 40xs/0.65, 100xs/1.25(Oil)					
Stage	Double Layer Mechanical Stage, 140*140mm, Moving Range 50*75mm					
Focusing	Low Position Coaxial Coarse & Fine Focusing. Fine Focusing Scale 0.002mm, With Tension Adjustment, With Up Limit					
Condenser	N.A.1.25 Abbe Condenser, Iris Diaphragm, With Filter Holder					
Filter	Blue, Yellow, Green					
Light Source	Plan-Concave Mirror					
	Built-in Illumination 6V20W Halogen Lamp					
A11.107 Biological Microscope Optional Accessories						
Dark Field	Dark Field Condenser, N.A. 0.7~0.9, Dry					A5D.1006-A
	Dark Field Condenser, N.A. 1.25~1.35, Immersion					A5D.1006-B

- 195mm Achromatic Objectives, Mechanical Tube Length 160mm
- WF10x/18mm, WF16x/11mm Wide Field Eyepiece
- Monocular, Binocular, Trinocular Head Optional
- Low Position Coaxial Coarse & Fine Focus System
- Double Layer Mechanical Stage




A11.1007-17 Biological Microscope



Sliding Binocular Head, Diopter Adjustable



Quadruple Nosepiece, Achromatic Objectives



Built-in Illumination, Halogen Lamp 6V/ 20W



A11.1007-16



A11.1007-27




A11.1007-17



A11.1007-17E



A11.1007-17C



A5F.1001 Epi-Fluorescent Kit



A11.1007 + A5F.1001 Epi-Fluorescent Kit



Opto-Edu (Beijing) Co., Ltd.



0086 13911110627



sale@optoedu.com



cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China