



A11.1009 1000X LED Light Source Biological Microscope For Student / School

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A11.1009
- Minimum Order Quantity: 1 pc
- Price: Negotiation
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Applications: Student|school
- Certification: CE|Rohs
- Drawtube: Binocular|Monocular
- Illumination: LED
- Product Name: Microscope
- Theory: Biological Microscope
- Highlight: **high powered microscope, binocular biological microscope, biological compound microscope**



Product Description

A11.1009 1000X LED light source Biological Microscope

High Quality Compact Design Student Biological Microscope;
WF10x/18mm Wide Field Eyepiece;
Monocular, Binocular Head 45° Inclined, 360° Rotatable;
Double Layer Mechanical Stage 110x120mm Moving 70x20mm;
Coaxial Coarse & Fine Focusing Konb, Light Source Brightness Adjustable.



OPTO-EDU (BEIJING) CO., LTD.

OPTO-EDU

F-1501 Wanda Plaza, No.18 Shijingshan Road, Beijing 100043, China
Tel:+8610 88696020 Fax:+8610 88696085


A11.1009

Student Biological Microscope





A11.1009-D Monocular Head



A11.1009-E Binocular Head




Achromatic Objectives



Double Layer Mechanical Stage



Brightness Adjustable



Abbe Condenser N.A.1.2, With Iris Diaphragm



A11.1009-D



A11.1009-E

Items	A11.1009 Student Biological Microscope	M	B
Head	Monocular Head, 45° Inclined, 360° Rotatable	●	
	Sliding Binocular Head, 45° Inclined, 360° Rotatable		●
Eyepiece	Wide Field WF10x/18mm	●	●
	WF5x	○	○
	WF16x	○	○
	WF10x Eyepiece With Pointer	○	○
Nosepiece	Triple Nosepiece	●	
	Quadruple Nosepiece		●
Achromatic Objective	4x/0.10, 10x/0.25, 40x/0.65(S)	●	
	4x/0.10, 10x/0.25, 40x/0.65(S), 100x/1.25(S,Oil)		●
	20x 60x	○	○
Magnification	40x~400x	●	
	40x~1000x		●
Stage	Double Layer Mechanical Stage, 110*120mm, Moving Range 70*20mm	●	●
Focusing	Coaxial Coarse & Fine Focusing	●	●
Condenser	Abbe Condenser N.A.1.2, With Iris Diaphragm, With Blue Filter, Rack & Pinion Up/Down Adjustment	●	●
Light Source	Incandescent Lamp Illumination 15W, Brightness Adjustable	●	
	LED Lamp Illumination 1W, Brightness Adjustable		●
	Mirror	○	○
Spare Parts	Dust Cover, Spare Bulb & Fuse	●	●

Note: "●" In Table Is Standard Outfits, "○" Is Optional Accessories.



Opto-Edu (Beijing) Co., Ltd.



0086 13911110627



sale@optoedu.com



cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China