



High Precision Monocular Compound Microscope Clip Stage WF10x Eyepiece A11.1141

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A11.1141
- Minimum Order Quantity: 1 pc
- Price: Negotiation
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Eyepiece: WF10x
- Objective: Achromatic 4x, 10x, 40x
- Nosepiece: Backward Triple Revolving Nosepiece
- Stage: Clip Stage
- Focusing: Coaxial Coarse & Fine Focusing Knob
- Illumination: Up & Bottom LED Rechargeable
- Highlight: **high powered microscope ,
binocular biological microscope ,
biological compound microscope**

Product Description

Achromatic 4x,10x,40x A11.1141 Up & Bottom LED Rechargeable Biological Microscope OPTO-EDU

Specifications:

A11.1141 Student Biological Microscope		M	B
Head	Monocular Head	o	-
	Binocular Head	-	o
Eyeiece	WF10x	o	o
Objective	Achromatic 4x,10x,40x	o	o
Nosepiece	Backward Triple Revolving Nosepiece	o	o
Stage	Clip Stage	o	o
Focusing	Coaxial Coarse & Fine Focusing Knob	o	o
Illumination	Up & Bottom LED Rechargeable	o	o





Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China