



## 2mm - 60mm View Field Bullet Comparison Microscope Digital Camera A18.1827

### Our Product Introduction

#### Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A18.1827
- Minimum Order Quantity: 1 pc
- Price: Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



#### Product Specification

- Total Magnification: 3.2x-192x
- Head: Binocular 45° Inclined, 48-75mm
- Eyepiece: WF10X/22mm, WF20X/12mm
- Step Zoom Lens: 0.8x, 1.25x, 2x, 3.2x, 4.8x
- View Field: Full Left, Full Right, Split Image, Overlapping Image, 2mm-60mm
- Auxiliary Lens: 0.4x, 2.0x
- Highlight: **bullet comparison microscope, binocular optical microscope**



Product Description

OPTO-EDU A18.1827 3.2~192x 3.0M 2mm-60mm View Field Digital Camera bullet comparison microscope

A18.1827 Comparsion Microscope	
T o t a l M a g n i f i c a t i o n	3.2x-192x
H e a d	Binocular 0°-45° Adjustable, 48-75mm
E y e p i e c e	WF10X/22mm, WF20X/12mm
S t e p Z o o m L e n s	0.8x, 1.25x, 2x, 3.2x, 4.8x
A u x i l i a r y L e n s	0.4x, 2.0x
V i e w F i e l d	Full Left, Full Right, Split Image, Overlapping Image, 2mm-60mm
W o r k i n g S t a g e	Moving Range: X-54mm, Y-54mm, Z-54mm, 0°-45° Inclination Adjustable Joint Left/Right Moving Range 62mm, Joint Up/Down Range 62mm Round Stage, 360° Rotating Bullet Holder 2 Sets
L i g h t S o u r c e	5W LED+12V/50W Halogen Coaxial Light Polarizing Light UV Filter Filter Red/Blue/Yellow/Green



Our Product Introduction

for more products please visit us on [cnoec.com](http://cnoec.com)

C a m e r a	3.0M Digital Camera
A d a p t e r	0.5x C-Mount
	Photo Adapter For Cannon SLR
F u n c t i o n P a n e l	Touch Type, User-Friendly Design, Easy To Operate And Accurate
P o w e r	Power Supply 80-240V

**OPTO-EDU**

OPTO-EDU (BEIJING) CO.,LTD.



**A18.1827**

**Comparsion Microscope**

#### HD Image

Visual HD images, static data acquisition via camera installation

#### Application

It is the main scientific and technical means of national criminal investigation technology, procuratorial technology and forensic identification, and has been widely used in the field of judicial identification technology.



0.4x Auxiliary Lens



Bullet Holder



Filter Red/Blue/Yellow/Green/UV



Polarizing Light Lens Set



Coaxial Light Attachment



WF10x & WF20x

### Coin Test Switching Field Of View



Split Image



Overlapping Image



**Opto-Edu (Beijing) Co., Ltd.**



0086 13911110627



sale@optoedu.com



cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China