



Digital Comparison Microscope

Our Product Introduction

for more products please visit us on cnoec.com

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A18.1831
- Minimum Order Quantity: 1 pc
- Price: Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Total Magnification: 3.2x-288x
- Head: Binocular 45° Inclined, 48-75mm
- Eyepiece: WF10X/22mm, WF20X/12mm
- Step Zoom Lens: 0.8x, 1.25x, 2x, 3.2x, 4.8x
- View Field: Full Left, Full Right, Split Image, Overlapping Image, 2mm-60mm
- Auxiliary Lens: 0.4x, 2.0x
- Highlight: **bullet comparison microscope, binocular optical microscope**



Product Description

OPTO-EDU A18.1831 3.2~288x Binocular Step Zoom forensic microscope

A18.1831 Comparsion Microscope

T
o
t
a
l
M
a
g
n
i
f
i
c
a
t
i
o
n

3.2x-288x

H
e
a
d

Binocular 45° Inclinded, 48-75mm

E
y
e
p
i
e
c
e

WF10X/22mm, WF20X/12mm

S
t
e
p
Z
o
o
m
L
e
n
s

0.8x, 1.25x, 2x, 3.2x, 4.8x

A
u
x
i
l
i
a
r
y
L
e
n
s

0.4x, 2.0x

V
i
e
w
F
i
e
l
d

Full Left, Full Right, Split Image, Overlapping Image, 2mm-60mm

W
o
r
k
i
n
g
S
t
a
g
e

Moving Range: X-55mm, Y-55mm, Z-65mm, 0°-45° Inclination Adjustable

Joint Left/Right Moving Range 62mm, Joint Up/Down Range 80mm

Round Stage, 360° Rotating

Bullet Holder 2 Sets

B
r
i
d
g
e

Motorized Up/Down

Light Source	5W LED+12V/50W Halogen+108 LEDS Ring Light UV Filter Filter Red/Blue/Yellow/Green
Camera	3.0M Digital Camera, Built-in CCD
Adapter	0.5x C-Mount
Photo Adapter	Photo Adapter For Cannon SLR
Function Panel	Touch Type, User-Friendly Design, Easy To Operate And Accurate
Power	Power Supply 80-240V

OPTO-EDU

OPTO-EDU (BEIJING) CO.,LTD.



A18.1831

Comparsion Microscope

HD Image

Visual HD images, static data acquisition via camera installation

Application

It is the main scientific and technical means of national criminal investigation technology, procuratorial technology and forensic identification, and has been widely used in the field of judicial identification technology.



WF10X/22mm, WF20X/12mm



Bullet Holder



Filter Red/Blue/Yellow/Green/UV

Coin Test Switching Field Of View



Split Image



Overlapping Image



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China