



A18.1822 5 Step Zoom Lens Forensic Comparison Microscope Motorized 3.2x - 320x Magnification

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A18.1822
- Minimum Order Quantity: 1 pc
- Price: Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Total Magnification: 3.2x-320x
- Eyepiece: WF10X/22mm, WF16X/13mm, WF20X/12mm
- Zoom: Dual Bridge Zoom 1x/1.6x
- View Field: Full Left, Full Right, Split Image, Overlapping Image, 2mm-60mm
- Head: Binocular 45° Inclined, 48-75mm
- Auxiliary Lens: 0.4x.2.0x
- Bridge: Motorized Joint Up/Down Range 0-150mm
- Highlight: **bullet comparison microscope, binocular optical microscope**



Product Description

OPTO-EDU A18.1822 5 Step Zoom Lens motorized comparison microscope

A18.1822 Motorized Comparison Microscope	
T o t a l M a g n i f i c a t i o n	3.2x-320x
H e a d	Binocular 45° Inclinded, 48-75mm
E y e p i e c e	WF10X/22mm, WF16X/13mm, WF20X/12mm
Z o o m	Dual Bridge Zoom 1x/1.6x
O b j e c t i v e	When The Bridge Zoom Is 1x, 5 Step Zoom Lens 0.8x,1.25x,2x,3.2x,4.8x
	When The Bridge Zoom Is 1.6x, 5 Step Zoom Lens 1.28x,2x,3.2x,5x,8x
A u x i l i a r y L e n s	0.4x.2.0x
V i e w F i e l d	Full Left, Full Right, Split Image, Overlapping Image
W o r k i n g S t a g e	Manual Adjustment Range: X-55mm, Y-55mm, Z-55mm, Round Stage Size:Φ90mm
	Motorized Adjustment Range: X-50mm, Y-50mm, Z-50mm
	Joint Left/Right Moving Range 54mm, Joint Up/Down Range 60mm
W o r k i n g D i s t a n c e	102mm

B r i d g e	Motorized Joint Up/Down Range 0-150mm
L i g h t S o u r c e	Dual Light 5W LED Coaxial Light Polarizing Light UV Filter And Red/Yellow/Blue/Green Filter
C a m e r a	3.0M Digital Camera, Photo Adapter For Cannon SLR
A d a p t e r	0.5x C-Mount
F u n c t i o n P a n e l	Touch Type, User-Friendly Design, Easy To Operate And Accurate
P o w e r	Power Supply 80-240V

OPTO-EDU

OPTO-EDU (BEIJING) CO.,LTD.



A18.1822

Motorized Comparison Microscope

Imaging System

A unique three-stage amplification for more accurate, faster and wider zooming

Application

It is the main scientific and technical means of national criminal investigation technology, procuratorial technology and forensic identification, and has been widely used in the field of judicial identification technology.



0.4x Auxiliary Lens



Polarizing Light Lens Set



Coaxial Light Attachment

Product Real Shot

■ Coin Test Switching Field Of View



Split Image



Overlapping Image



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China