



Desktop Fingerprint Forensic Comparison Microscope A18.1839 Focusing Range 90mm

Our Product Introduction


Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A18.1839
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Total Magnification: 6.5x/10x
- Screen Size: 192x165mm
- Focusing Range: 90mm
- View Field: 29.5x25.5/19.4x166mm
- Working Distance: 80/47mm
- Light Source: 220V/50HZ, Halogen Lamp 300W/120V
- Highlight: **bullet comparison microscope, binocular optical microscope**



Product Description

OPTO-EDU A18.1839 Screen Size 192x165mm Focusing Range 90mm Desktop Fingerprint Comparator

A18.1839 Desktop Fingerprint Comparator	
Total Magnification	6.5x/10x
Screen Size	192x165mm
Focusing Range	90mm
View Field	29.5x25.5/19.4x166mm
Working Distance	80/47mm
Light Source	220V/50HZ, Halogen Lamp 300W/120V



Opto-Edu (Beijing) Co., Ltd.



0086 13911110627



sale@optoedu.com



cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China