



Multi Functional Micro Trace Forensic Comparison Microscope A18.1840 0.7 - 4.5x

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A18.1840
- Minimum Order Quantity: 1 pc
- Price: Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Total Magnification: 0.7~4.5x
- View Field: 72mm
- Manual Focus: 70mm
- Stand: Pole Stand Up/Down 200mm
- Camera: 3.0M Digital Camera
- Light Source: UV LED5W
- Other Accessories: Tripod Stand, Computer, Incl. Laptop
- Highlight: **bullet comparison microscope, binocular optical microscope**



Product Description

OPTO-EDU A18.1840 0.7~4.5x View Field 72mm Multifunctional Micro Trace Identification Instrument

A18.1840 Multifunctional Micro Trace Identification Instrument	
Total Magnification	0.7~4.5x
View Field	72mm
Manual Focus	70mm
Stand	Pole Stand Up/Down 200mm
Camera	3.0M Digital Camera
Light Source	UV LED5W
Other Accessories	Tripod Stand, Computer, Incl. Laptop



Opto-Edu (Beijing) Co., Ltd.



0086 13911110627



sale@optoedu.com



cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China