



A18.1841 12V 20W Microscope Halogen Lamp Brightness Adjustable Side Halogen Light

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A18.1841
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Light Source: 12V20W Halogen Lamp
- Colour: Black
- Use: Microscope Assisted Light Source
- Highlight: **bullet comparison microscope, binocular optical microscope**



for more products please visit us on cnoec.com

Product Description

OPTO-EDU A18.1841 12V20W Halogen Lamp Brightness Adjustable Side Halogen Light

Side light

Voltage input: 220v

Light source: Halogen bulb 12v20w, adjustable brightness, hose design on the neck, 360° rotation adjustable

Main scope of application: as an auxiliary light source for indentation detection, visual illumination source of the machine, etc.



Opto-Edu (Beijing) Co., Ltd.



0086 13911110627



sale@optoedu.com



cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China