



## Binocular Stereo Microscope 20X 40X

### Our Product Introduction

for more products please visit us on [cnoec.com](http://cnoec.com)

#### Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE,
- Model Number: A22.1233
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



#### Product Specification

- Head: Rotatable 360°, Inclination 45° Binocular, With Locked-in Eyepiece, Left & Right Diopter Adjustable  $\pm 5^\circ$ , Head Interpupillary 55-75mm
- Eyepiece: Extra Wide Field, Plan WF10x/22mm, High Eyepoint, Eyepiece Tube  $\Phi$  30mmWF10X
- Objective: Three Kind Of Objective 90° Rotatable For Selection: 1X2X, 1X3X, 2X4X
- Working Stage:  $\Phi$ 60 White/Black And  $\Phi$ 60 Frosted Glass
- Power Supply: AC 100V-240V Wide Range, CE Approval
- Light Source: No Light Source Or Reflect 1WLED, Transmit 7x20mA
- Highlight: **stereo zoom microscopes,**  
**stereo binocular microscopes,**  
**stereo dissecting microscope**



## Product Description

### OPTO-EDU Cheap Student Boom Stand Inspection Electron Pcb Repair Optical 20x 40x Binocular Stereo Microscope Sunshine

Stereo microscope is widely used in teaching demonstration, bio-engineering, IT industry testing and other fields. A22.1233 has a high resolution and large field of view clarity, so it is also good at electronics industry and precision machinery industry parts assembly and testing, agricultural seed inspection and so on.

A22.1233-PA Pole Stand+Touch Switch



A22.1233-PB Pole Stand+Dial Switch



A22.1233-TA Track Stand + Touch Switch



A22.1233-TB Track Stand+Dial Switch



Model	P	PA	PB	TA	TB
Head	Head Fixed, Inclination 45° Binocular, With Locked-in Eyepiece, Left & Right Diopter Adjustable ±5°, Head Interpupillary 55-75mm				
Eyepiece	WF10X				
Objective	Three Kind Of 180°Rotatable Objective For Selection: 1X2X, 1X3X, 2X4X				
Light Source	No Light Source	Reflect 1WLED, Transmit 7x20mA			
Stand	Pole Stand	Pole Stand Touch Switch	Pole Stand Dial Switch	Track Stand Touch Switch	Track Stand Dial Switch
Working Stage	Φ60 White/Black	Φ60 Frosted Glass			
Power Supply	AC 100V-240V Wide Range, CE Approval				
Optional Accessories					
Eyepiece	WF5X				
	SWF10X				
	WF15X				
	WF20X				

Working Stage	Φ60 White/Black
---------------	-----------------



**Opto-Edu (Beijing) Co., Ltd.**



0086 13911110627



sale@optoedu.com



cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China