



Cnoec A21.3601-B3 4.5x Objective Stereo Optical Microscope

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE,
- Model Number: A21.3601-B3
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Certification: CE
- Drawtube: Monocular
- C-Mount: 0.5x, W.D.100mm
- Objective: 0.7x-4.5x
- Magnification: 0.35x-2.25x, Zoom Power
- Focusing Type: B3 Pillar Stand, Column Diameter $\Phi 25\text{mm}$, Height 323mm Or 270mm
- Stand: B3 Pillar Stand, Column Diameter $\Phi 25\text{mm}$, Height 323mm Or 270mm
- Highlight: **4.5x Stereo Optical Microscope,**
4.5x video microscope,
cnoec Stereo Optical Microscope



Product Description

Stereo Optical Microscope Monocular Zoom Video Microscope A21.3601-B3

Items	A21.3601 Mono Zoom Video Microscope	A21.3601	A21.3601-B3	A21.3601-STL7
Head	Monocular	•	•	•
C-Mount	0.5x, W.D.100mm	•	•	•
Objective	0.7x-4.5x	•	•	•
Magnification	0.35x-2.25x, Zoom Power	•	•	•
Focusing Type	XLB-A1 Basket Holder, Range 62mm, Mirror Aperture Φ50mm	•	•	•
Stand	None, Main Body Only	•		
	B3 Pillar Stand, Column diameter Φ25mm, Height 323mm or 270mm		•	
	STL7 Stand, Diameter Φ25mm, Height 340mm, Size 380x260x20mm			•
Items	A21.3601 Mono Zoom Video Microscope Optional Accessories			
C-Mount	0.35x			
	1.0x			
Auxiliary Objective	0.5x			
	2.0x			
Light Source	LED Ring Lamp			

Eyepiece	Standard Objective 0.7x-4.5x		Auxiliary objective 0.5x		Auxiliary objective 2.0x	
	W.D.100mm		W.D.198mm		W.D.46mm	
	Magnification	Field Of View	Magnification	Field Of View	Magnification	Field Of View
0.35x	0.245x	23.0	0.1225x	46.0	0.49x	11.6
	1.575x	3.0	0.7875x	7.0	3.15x	1.7
0.5x	0.35x	17.6	0.175x	35.0	0.7x	8.8
	2.25x	2.7	1.125x	5.0	4.5x	1.3
1.0x	0.7x	9.0	0.35x	17.6	1.4x	4.0
	4.5x	1.6	2.25x	3.0	9x	0.7



OPTO-EDU (BEIJING) CO., LTD.

OPTO-EDU

F-1501 Wanda Plaza, No.18 Shijingshan Road, Beijing 100043, China
Tel:+8610 88696020 Fax:+8610 88696085

A21.3601 Mono Zoom Video Microscope



Product Overview



Product Specification



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China