



Lab OPTO-EDU A11.2601-T Trinocular Biological Microscope

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A11.2601-T
- Minimum Order Quantity: 1 pc
- Price: Negotiation
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Head: Seidentopf Trinocular Head Inclined 30°
- Applications: Student|school
- Certification: CE|Rohs
- Drawtube: Trinocular
- Illumination: Halogen Lamp|LED
- Max Magnification: 1000X
- Product Name: Microscope
- Highlight: **Trinocular Biological Microscope, OPTO-EDU lab microscope, A11.2601-T Biological Microscope**



More Images



Product Description

A11.2601-T Student Microscope.

[illegible]



Biological Microscope CNOPTEC Series

A11.2601

- ♦ Wide Field Eyepiece WF10x/18mm
- ♦ High Contrast Achromatic Objective 4x, 10x, 40x, 100x
- ♦ Mechanical Stage 115x125mm, Moving Range 70x35mm
- ♦ Full Metal Body High Rigid Construction
- ♦ 3W LED, Cordless Rechargeable Function Optional



Pretty

Durable

Good Effect

A11.2601 Biological Microscope Is Designed For Teaching, Clinical Verify. The Series Microscope With Original Style, Steady Structure, Convenient Operation And Clear Image Is Suitable For Observing Various Biological Specimens, And Mostly Applied In School, College, Hospital And Family.



01 Head

WF10x-18mm
Eyepiece Is Anti-fungal, Anti-fungal System
Built In Head

02 Objective

Achromatic Objective: 4x, 10x, 40x(s),
100x/1.25(oil)(s)
Objective Is Anti-fungal



03 Stage

Mechanical Stage
Stage Size:125mm×115mm,
Travel: 35mm×75mm

04 Focus

Coaxial Coarse & Fine Focus
Stage Moving Knob Can Move Stage
Turn Left Or Right



05 Condenser

Abbe Condenser
N.A. 1.25, Adjustable Aperture, Aperture
Center Is Adjustable
Built In Blue Filter, $\phi 38.5\text{mm}$

06 Light Source

LED Illumination Systems (Non-recharge-
able)
Rotation Potentiometer With Power Switch





Opto-Edu (Beijing) Co., Ltd.



0086 13911110627



sale@optoedu.com



cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China