



12.5X Medical 30mm Dental Lab Microscope Surgical Operating

Our Product Introduction

for more products please visit us on cnoec.com

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE,
- Model Number: A41.1903
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month





Product Specification

- Binocular Observation: Tilting Head 0°-180° Viewing Angle Adjustable For User
- Magnification Of Binocular: 8X
- Magnification Of Eyepiece: 12.5X
- Interpupillary Distance: 50mm-80mm
- Diopter: ±7D
- Magnification Changer: 3-step Magnification Change: 0.6X, 1X, 1.6X
- Focal Length Of Objective: Two Objectives: F=200 & F=250mm (M45x0.75mm)
- Total Magnification: 4.8X, 8X, 12.8X, 6X, 10X, 16X
- Linear Field: 50mm, 30mm, 19mm, 40mm, 24mm, 15mm
- Filter: Built-in Green And Yellow Filters
- Highlight: **12.5X dental lab microscope, 30mm dental lab microscope, 30mm Surgical Operating microscope**



More Images



Product Description

A41.1903 Series Operation Microscope Specification

A41.1903 Series Operating Microscopes, is the basic economical high quality series. It use fibre-optic cable for coaxial illumination, with 3 steps zoom body, has total magnification 2.3x~12x, featured with flexible handling and easy operation. It is widely used in hospital for ENT (ear-nose-throat), dental, ophthalmology, gynecology, surgery, orthopedics and etc.

OPTO-EDU

A41.1903 Operating Microscope

Dental, ENT, Ophthalmology
Gynecology



Tilting Head 0° - 180°
Viewing Angle Adjustable



Digital Camera & Adapter

Demonstrator (Optional)

Full Range Accessories Optional

Camera Adapter | Beam Splitter 50:50 | Objective

 <p>WF10x/23mm Eyepiece</p>	 <p>SLR Camera Adapter</p>	 <p>LWD Objective 300mm</p>
 <p>Teaching / Assistant Head</p>	 <p>50:50 Beam Splitter</p>	 <p>Digital Camera Adapter</p>



Model	A41.1903
Binocular Head	Tilting 0~180°
Main Usage	Dental, ENT, Ophthalmology, Gynecology
Total Magnification	With Objective F=200mm: 6.0x, 10x, 16x With Objective F=250mm: 4.8x, 8.0x, 12.8x
Total Magnification (Optional)	With Objective F=300mm: 4.0x, 6.7x, 10.7x With Objective F=400mm: 3.0x, 5.0x, 8.0x
Eyepiece	WF12.5x/20mm
Head Magnification	8x
Interpupillary Distance	50~80mm Adjustable
Diopter	Both Left & Right Eyetube Diopter Adjustable $\pm 7D$
Zoom Body	3 Steps Zoom: 0.6x, 1.0x, 1.6x
Objective	F=200mm For Ear, Nose, Ophthalmology F=250mm For Dental
Objective (Optional)	F=300mm For Ear, Nose F=400mm For Throat

Focusing	Fine Focusing Range 10mm	
Filter	Green, Yellow Filters Built-in	
Balancing Arm	2-Part Arm With Universal Joints, Counterweight Adjustable And Lockable	
Stand	2-Part Column + Five Star Base With Caster	
Illumination System	Coaxial Illumination, 10W LED, Brightness Adjustable, >40000 Lux	
Power Supply	Wide Range AC100~240V, DC12V 3A, 36W	
A41.1903 Series Operation Microscope Optional Accessories		
Eyepiece	WF10x/23mm	A51.1901-1023
Objective	F=300mm	A52.1901-300
	F=400mm	A52.1901-400
	F=200mm/250mm/300mm Zoom, Connection Thread: M47x0.5mm(Zeiss) / M65x0.5mm(Leica), Or M47x0.75mm / M51x0.75mm / M56x0.75mm, Or Customized According To Customer Requirements	A52.1902-01
	F=200mm/250mm/300mm Zoom, With Splash Cover Connection Thread: M47x0.5mm(Zeiss) / M65x0.5mm(Leica), Or M47x0.75mm / M51x0.75mm / M56x0.75mm, Or Customized According To Customer Requirements	A52.1902-02
	50:50 Beamsplitter	A55.1906
Adapter	C-Mount 0.3x For 1/3" Digital Camera	A55.1903
	SLR Camera Adapter For NIKON DSLR Or CANON SLR	A55.1904
	Teaching Head, Demonstrator	A55.1910
Head	Tilttable Head, Inclinable Angle: 0-180°, Eyepiece 12.5X/20mm, Interpupillary Distance 50~80mm, Magnification 8x, F=160mm, Distance of Center 22mm. Optional Eyepiece: 10x/23mm	A53.1903-01
	Tilttable Head, With PD Adjust Knob Inclinable Angle: 0-180°, Eyepiece 12.5X/20mm, Interpupillary Distance 50~80mm, Magnification 8x, F=160mm, Distance of Center 22mm. Optional Eyepiece: 10x/23mm	A53.1903-02
Mount	Table Mount Clamp	A54.1901-T
	Wall Mount Bracket	A54.1901-W
Monitor Mount	Monitor Mount: 1 Fix Arm And 1 Gas Spring Arm With Tension Adjustable. Band Clamp. Fit 15-27" Screen. VESA Standard 75x75xmm & 100x100mm. Load Capacity 0~8 Kgs, Swivel 180°, Tilt -80°~+90°, Rotation 360°	A54.1910-112B
	Monitor Mount: 1 Fix Arm And 1 Mechanic Spring Arm With Tension Adjustable. Band Clamp. Fit 15- 27" Screen. VESA Standard 75x75xmm & 100x100mm. Load Capacity 0~8 Kgs, Swivel 180°, Tilt -80°~+90°, Rotation 360°	A54.1910-312B



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China