



1000X 5.0M Resolution 9.7" Portable Lcd Digital Microscope

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A33.1502
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Authenticate: CE
- Export Area: Europe & North America
- Monthly Sales: 5000 Units
- Application Field: Efficient And Laboratory
- Magnification: 1000x
- Resolution: 5.0M
- Highlight: 9.7" Portable Lcd Digital Microscope,
1000X Portable Lcd Digital Microscope,
5.0M Portable Lcd Digital Microscope



More Images



Product Description

1. 9.7" LCD, 5.0M, Trinocular, WF10x/20, Quaduple, Achromatic HC 4x,10x,40x100x, Double Layer Mechanical Graphite 143*132mm, Moving 75x50mm, Abbe N.A.1.25, Coaxial, 3W LED, 1.0x C Mount
2. 1.9.7" LCD Touch Screen Pad, With Android System + Measuring Software
3. 2.High Eyepoint Wide Field Eyepiece WF10x/20mm Dia.30mm
4. 3.Achromatic HC 4x,10x,40x100x Objective, Upgradeable To Infinity
5. 4.Double Layer Mechanical Stage 143*132mm, Upgradeable To 216x150mm
6. 5.3W LED Kohler Illumination With Removable Light House For Easy Maintenance

OPTO-EDU

A33.1502 Research Level Digital Microscope

9.7" LCD | Android Pad | 5.0M



9.7" Android Pad + 5.0M Camera


Sharp Image | All-in-One | Photo & Video



Explore Micro World Together

Multipul Usage | Easy To Share | Nice For Teaching





A33.1502 9.7" LCD Digital Biological Microscope	
Head	Trinocular Head, 30° Inclined, Interpupillary Distance 48-75mm, Both Eyepiece Tube Diopter Adjustable
Eyepiece	WF10x/20mm, Diopter Adjustable, Dia. 30mm
Objective	High Contrast Chromatic Free 4x, 10x, 40x(S), 100x(Oil, S)
Focus Knob	Coaxial Coarse & Fine Focusing, Fine Division 0.002mm, Coarse Stroke 37.7mm/Rotation, Focusing Range 24mm
Nosepiece	Quadruple, Backward With Click Stops
Stage	Double Layer Mechanical Stage 142x132mm, Graphite Surface, Moving Range 75x50mm
Condenser	Abbe Condenser N.A.1.25, Center Adjustable
Light Source	3W LED Brightness Adjustable
Adapter	C-Mount 1.0x
Package	53*32.5*53.5cm/11.5Kg, 1 Set/Carton
9.7" LCD Pad Digital Camera (A59.3520)	
LCD	9.7" HD IPS Touch Screen With Backlight
Camera	5.0M 1/2.5" CMOS Digital Camera
Resolution	2048x1536 @15FPS
Pixel	2.2umx2.2um
Data Storage	TF Card or Built-in Memory
Data Output	USB2.0, HDMI, TF Card
Photo	5.0M .JPEG
Video	1080p Video Record
System	Android 5.1 System
	90 Languages
	Support Mouse + Keyboard Control
	Support WIFI
Function	Support Bluetooth
	S-Eye Software & Instruction Manual Installed
	Mark: Point, Crosshair, Coordinate, Text
	Length Measure: Straight Line Length, Broken Line Length, Curve Length, Parallel Line Distance, Point To Line Distance
	Distance Between 2 Circle Centers: Radius Fixed Circle, Two-point Fixed Circle, Three-point Fixed Circle
	Geometric Measure: Line Segment Length, Radius Fixed Circle, Two-point Fixed Circle, Three-point Fixed Circle, Concentric Circle
	Geometric Area Measure: Polygon, Square
Power	DC 12V/2A Power Adaptor, AC 100~240V
Size	Pad 238x182x12mm, Camera 40x39mm





0086 13911110627



sale@optoedu.com



cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China