



## OPTO-EDU 7x-45x Electron Stereo Optical Microscope

Our Product Introduction

for more products please visit us on [cnoec.com](http://cnoec.com)

### Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE,
- Model Number: A23.3645N-R2H1
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month




### Product Specification

- LCD Camera: 9" LCD Screen + Digital Camera
- Head: Trinocular Head, Light Split Fixed 50:50
- Eyepiece: WF10x/20mm
- Zoom Ratio: 1:6.5
- Zoom Lens: 0.7~4.5x
- Magnification: 7x-45x
- Working Distance: 100mm, With Auxiliary Objectives 30~165mm
- Field Of View: 1~95.2mm
- Stand: Pole Stand, Size 200\*255\*22mm
- Highlight: Electron Stereo Optical Microscope, 7x Stereo Optical Microscope, OPTO-EDU Microscope



### More Images



## Product Description


## OPTO-EDU A23.3645N-R2H1 7-45x Zoom Stereo Microscope

A23.3645N-R2 Series DIY Phone PCB Repair Microscope		A23.3645N-R2	A23.3645N-R2L1	A23.3645N-R2H1
Head	Trinocular Head, With Eyepiece Adapter	●		
	LCD Screen 9" With Digital Camera, AV Output		●	
	HDMI+USB Digital Camera 38.0M Pixels			●
	45° Inclined, 360°Rotatable, Both Tube ±5 Diopter Adjustable, Interpupillary Adjustment 54-75mm			
Eyepiece	WF10x/20mm, Both With Eye-cup			
Zoom Lens	0.7~4.5x			
Zoom Ratio	1:6.5			
Magnification	7x-45x, Up To 3.5x-180x With Optional Eyepieces And Auxiliary Lens			
W.D.	Working Distance100mm, With Standard 1x Lens			
F.O.V	Field Of View 1~95.2mm			
Stand	Large Pole Stand, Size 200*255*8mm, Pole Height 240mm, Dia.32mm, With Heat Resist Soft Rubber Pad, Working Area For Phone PCB Repair			
Light Source	Upper 56 LEDs Right Light Brightness Adjustable			
Optional Accessories			Cata.No.	
Eyepiece	WF15x/15mm, High-eyepoint	A51.3621-1515		
	WF20x/10mm, High-eyepoint	A51.3621-2010		
	WF25x/9mm, High-eyepoint	A51.3621-2509		
	WF30x/8mm, High-eyepoint	A51.3621-3008		
Auxiliary Lens	0.5X,W.D. 165mm	A52.3621-0.5x		
	1.5X,W.D. 45mm	A52.3621-1.5x		
	2.0X,W.D. 30mm	A52.3621-2.0x		

A23.3645N Zoom Stereo Microscope Optical Data Sheet									
Objective	Standard 1.0x		Auxiliary 0.5X		Auxiliary 1.5X		Auxiliary 2X		
W.D.	100mm		165mm		45mm		30mm		
Eyepiece	Mag	View Field	Mag	View Field	Mag	View Field	Mag	View Field	
10X/20mm	7.0X	28.6mm	3.5X	57.1mm	10.5X	19mm	14X	14.3mm	
	45.0X	4.4mm	22.5X	8.9mm	67.5X	3mm	90X	2.2mm	
15X/15mm	10.5X	21.1mm	5.25X	42.8mm	15.75X	14.3mm	21X	10.7mm	
	67.5X	3.3mm	33.75X	6.7mm	101.25X	2.2mm	135X	1.7mm	
20X/10mm	14.0X	14.3mm	7.0X	28.6mm	21.0X	9.5mm	28X	7.1mm	
	90.0X	2.2mm	45.0X	4.4mm	135.0X	1.5mm	180X	1.1mm	



Upgrade With HDMI  
Or LCD Digital  
Camera



Featuring Simul-Focal Head,  
View Image In Eyepiece &  
Camera At Same Time



56 LEDs Ring Light  
Brightness Adjustable

## Specially Designed For Phone Repair Work!



Soft Heat Resist Rubber Pad  
Multiple Storage Area For Tools



## All Kinds Models Optional



A23.3645N-R2



A23.3645N-R2H1



A23.3645N-R2L1



BACK HOME



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China