



4K HDMI USB3.0 Output Tf Card Camera 12M 60FPS

Our Product Introduction

for more products please visit us on cnoec.com

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A59.3513
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month

OPTO-EDU



Product Specification

- Sensor: IMX226
- Pixel: 12M
- Resolution: 3840*2160@60FPS
- Output: HDMI2.0+USB3.0 Sync Output
- Pixel Size: 1.85μm*1.85μm
- Power: DC-12V/1A
- Highlight: 4K tf card camera, HDMI tf card camera, USB3.0 tf card camera

OPTO-EDU



More Images

OPTO-EDU



Product Description

5.0M, 9" LCD, 1/2.8" CMOS, Screen Resolution 1024*600, USB2.0 Real Time Image Output, TF Card, USB Mouse Control Menu For Measure, DC12V/2A, 100~240V.
4K HDMI + USB 3.0 Image Output + TF Card Storage
12M, 3840x2160@60FPS CMOS Sony Sensor
USB3.0 Real Time Image Output + S-Eye Measure Software
USB Mouse Control Menu On HDMI Output For Measure Function
DC12V/2A, 100~240V Power Adapter



OPTO-EDU (BEIJING) CO., LTD.

OPTO-EDU


F-1501 Wanda Plaza, No.18 Shijingshan Road, Beijing 100043, China
Tel:+8610 88696020 Fax:+8610 88696085

A59.3513


4K HDMI + USB 3.0 Digital Camera, 12M




Product Details





HDMI Mouse Menu



Product Size



S-EYE Measure Software For Windows



Items	A59.3513 4K HDMI Industry Digital Camera
Sensor	IMX226
Pixel	12M
Resolution	3840*2160@60FPS
Output	HDMI2.0+USB3.0 sync output
Pixel size	1.85μm*1.85μm
Power	DC-12V/1A
Function	Basic Function: Image Capture, Video, Measurement, Image Compare
	Image Adjustment: Exposure, Color Temperature, White Balance
	Camera Function: Mirror, Image Freeze, Flip, Wide Dynamic, Zoom In, Zoom Out
	Grid 8 Sets of Vertical/Horizontal Lines, Color/Size Can Be Adjusted
Package	Resolution 3840*2160
	Size 75*75*38.8mm
	Weight 0.2kg
	Accessory DC-12V/1A adapter



Opto-Edu (Beijing) Co., Ltd.

0086 13911110627

sale@optoedu.com

cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China

