



Binocular Led Student 18mm Eyepiece Biological Microscope

Our Product Introduction

for more products please visit us on cnoec.com

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A11.0216
- Minimum Order Quantity: 1 pc
- Price: Negotiation
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Eyepiece: WF10X(Φ 18mm)
- Eyepiece Tube: Binocular(Inclination Of 30°) 360° Rotatable
- Focus System: Coaxial Coarse/Fine Focus System With Up Stop, Minimum Division Of Fine Focusing: 2 μ m
- Nosepiece: Quadruple(Backward Ball Bearing Inner Locating)
- Stage: Double Layer Mechanical (Size:125mmX115mm, Move Range: 70mmX30mm)
- Abbe Condenser: N.A.1.25, Rack & Pinion Adjustable
- Filter: Frosted Filter
- Light Source: 1W LED With Brightness Control, Optional Rechargeable Battery, External 5V Power Adapter Or Portable Source
- Highlight: **binocular biological microscope, 18mm Eyepiece Biological Microscope, Led opto edu microscope**



More Images



Product Description

Binocular, WF10x/18, Quadruple, Achromatic 4x10x40x100x, Mechanical Stage 115x125mm, Abbe NA1.25, Coaxial Coarse & Fine Focusing, 1W LED Low Voltage
 Binocular Head With Left Eyepiece Tube Diopter Adjustable
 Achromatic/Plan Objective On Quaduple With Slot For Polarizer
 Colorful Handle For Easy Carrying Optional
 1W LED Light Source Brightness Adjustable
 Low Voltage Power Supply By External 5V Adapter, Or USB Power Bank


THE KEY TO MICRO WORLD

OPTO-EDU





A11.0216 Student Biological Microscope

NEW 2020!



Binocular Head, Left
Tube Diopter
Adjustable

High Quality
Achromatic Objective
4x, 10x, 40x, 100x



1W LED + Low
Voltage Power Supply
By External 5V
Adapter, Or USB
Power Bank



Product Details



Binocular Head Inclination 30°, Rotatable 360°



Quadruple Nosepiece With Achromatic Objectives



Slot For Polarizing Analyzer (Optional)



Double Layer Mechanical Working Stage



Abbe Condenser N.A.1.25, Rack & Pinion Adjustable



1W LED, Brightness Adjustable




Handle For Easy Carrying Optional



Low Voltage Power Supply By External 5V Adapter, Or USB Power Bank

A11.0216 Biological Microscope	
Eyeiece	WF10X(Φ18mm)
Objective(Standard)	Achromatic 4X/0.10
	Achromatic 10X/0.25
	Achromatic 40X/0.65(Spring)
	Achromatic 100X/1.25(Spring, Oil)
	Plan Achromatic 4X/0.10
Objective(Optional)	Plan Achromatic 10X/0.25
	Plan Achromatic 40X/0.65(Spring)
	Plan Achromatic 100X/1.25(Spring, Oil)
Eyeiece Tube	Binocular(Inclination Of 30°),360° Rotatable

Focus System	Coaxial Coarse/Fine Focus System With Up Stop, Minimum Division Of Fine Focusing: 2μm
Nosepiece	Quadruple(Backward Ball Bearing Inner Locating)
Stage	Double Layer Mechanical (Size:125mmX115mm, Move Range: 70mmX30mm)
Abbe Condenser	N.A.1.25, Rack & Pinion Adjustable
Filter	Frosted Filter
Collector	For LED Lamp
Light Source	1W LED With Brightness Control, Optional Rechargeable Battery, External 5V Power Adapter Or Portable Source



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China