



4.3" Screen 8.0m Resolution Led Portable Lcd Digital Microscope

Our Product Introduction

for more products please visit us on cnoec.com

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A33.1511
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Magnification: Total Magnification 70x~280x On LCD Screen
- LCD & Camera: 4.3" LCD Screen 8.0M Resolution
- Objective: Achromatic 10x/0.25, 20x/0.4x, 40x/0.65(Spring)
- Nosepiece: Outward Triple With Clip Stop
- Working Stage: Plain Working Stage 95x95mm
- Focusing: Coaxial Coarse & Fine Focusing Knob
- Highlight: **4.3" portable lcd digital microscope, 8.0m portable lcd digital microscope, led portable lcd digital microscope**



More Images



Product Description

4.3" LCD, Triple, Achromatic 10x20x40x, Vernial Scale Stage 95x95/50x30mm, Disc Diaphragm, Up/Bottom LED, With Acce. Set, With Carry Box

A33.1511 4.3" LCD Digital Microscope, 2.0M	
Magnification	Total Magnification 70x~280x on LCD Screen
LCD & Camera	4.3" LCD Screen 8.0M Resolution, With Built-in HD Digital Camera 2.0M All-in-One Design LCD Menu Support Brightness/Contrast/Color Adjust, Switch Between 16:9/4:3 Modes, Viewing Area 98*56mm, Support 10 Languages.
Objective	Achromatic 10x/0.25, 20x/0.4x, 40x/0.65(Spring)
Nosepiece	Outward Triple With Clip Stop
Working Stage	Plain Working Stage 95x95mm, With Vernial Scale Ruler Moving Range 50x30mm, With Disc Diaphragm Condenser
Focusing	Coaxial Coarse & Fine Focusing Knob, Focusing Range 16mm
Illumination	Dual Light Up/Bottom LED, Brightness Adjustable
Power	Powered By USB Adapter AC100V-240V 50/60HZ, or By USB Port of Computer, USB Bank For Outdoor Work
Accessory	Full Set Accessory Kit,
Packing	With Carry Box, Packed In Master Carton 25*31*43cm/3Kgs

A33.1511

4.3" LCD Digital Microscope, 2.0M



4.3" HD LCD
2.0M Digital Camera



Triple Nosepiece
Achromatic 10x20x40x
Clear Imaging



Dual Light
Up/Bottom LED



Product Details



LCD screen 45° tiltable, ergonomic design to meet observations from different directions easily.




X/Y axes moving ruler can quickly & accurately move the specimen into the viewing field




Elegant smooth design, heavy metal base supply steady stand for firm view.



Vernier scale stage, with single lens condenser, with disc diaphragm & different color filters.



Coaxial coarse & fine focusing, to focus the image clear easily



AC100V-240V 50/60HZ, external power supply without battery



Research



Pet Hospital



School



Aquaculture



Beauty mites



Children's Science



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China