



## Compound Rohs 1000x Dual Led Digital Biological Microscope

Our Product Introduction

for more products please visit us on [cnoec.com](http://cnoec.com)

### Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A11.1123
- Minimum Order Quantity: 1 pc
- Price: Negotiation
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



### Product Specification

- Applications: Student|school|Hobby|Home
- Certification: CE|Rohs
- Drawtube: Binocular
- Illumination: LED
- Max Magnification: 1000X
- Product Name: Microscope
- Theory: Biological Microscope
- Highlight: **dual led digital biological microscope,  
1000x digital biological microscope,  
Rohs biological compound microscope**




## Product Description

### Portable Cordless Biological Microscope 1000X For School / Home A11.1123

#### Description:

We have various Student Biological Microscope for your choice for student educational purpose. Designed with high quality optical system, steady and heavy materials, we can supply a real optical microscope for all students at very good prices. Sold to world market in large quantity, these student microscopes has become most popular tools in primary, middle and high school.

#### Product Details




Dual LED Light, Up/Bottom Illumination, Brightness Adjustable



Both Eyepiece Tube Diopter Adjustable



Coaxial Coarse And Fine Focusing Knobs



Double Layers Mechanical Stage



WF10x 2 pcs, 1pc With Pointer

## Product Display









A11.1123 Side View



A11.1123 Back View

## Product Specification

A11.1123 Economical Binocular Biological Microscope	
Head	Sliding Binocular Head, 45°Inclined, 360°Rotatable
Eyeiece	WF10x 2 pcs, 1pc With Pointer
Nosepiece	Quadruple Revolving Nosepiece
Objective	Achromatic 4x,10x,40x(Spring),100x(Spring,Oil)
Working Stage	Double Layers Mechanical Stage 100x90mm, Moving Range 30x60mm
Condenser	Abbe Condenser N.A1.25 With Iris Diaphragm
Focusing	Coaxial Coarse And Fine Focusing Knobs
Illumination	Dual LED Light, Up/Bottom Illumination, Brightness Adjustable
Power Supply	By External AC/DC Adaptor By 3xAA Rechargeable Batteries (If Shipped By Air, Battery Will NOT Be Included)



Opto-Edu (Beijing) Co., Ltd.



0086 13911110627



sale@optoedu.com



cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China