



A23.1503 Zoom Stereo Microscope 0.7-4.5x Mobile Phone Repair Can Up To 3.5x-180x

Our Product Introduction

for more products please visit us on cnoec.com

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A23.1503
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Head: HD+USB Digital Camera 38.0M Pixels
- Eyepiece: WF10x/20mm, Both With Eye-cup
- Zoom Lens: 0.7~4.5x
- Zoom Ratio: 1:6.5
- Magnification: 7x-45x, Up To 3.5x-180x With Optional Eyepieces And Auxiliary Lens
- W.D.: Working Distance 100mm, With Standard 1x Lens
- F.O.V: Field Of View 1~95.2mm
- Stand: Large Pole Stand, Size 200*255*8mm
- Light Source: Upper 56 LEDs Right Light Brightness Adjustable
- Highlight: 0.7x zoom stereo microscope, 4.5x zoom stereo microscope, zoom stereomicroscope



More Images



Product Description

Head: Trinocular Head, Both Eye Tube Diopter Adjustable
Light Split Switch E100:P0/E0:P100
With Eyepiece Adapter + 1.0x C-Mount

Eyepiece: WF10x/20mm, Both With Eye-cup

Zoom Lens: 0.7~4.5x


Zoom Ratio: 1:6.5

Magnification: 7x-45x, Up To 3.5x-180x With Optional Eyepieces And Auxiliary Lens


W.D.: Working Distance 100mm, With Standard 1.0x Lens


F.O.V: Field Of View 1~95.2mm

Stand: Large Pole Stand, Size 480x320mm, Pole Height 240mm, Dia.32mm, With Heat Resist Soft Rubber Pad, Working Area For Phone PCB Repair



Microscope For Mobile Phone Repair





Large Size Working Area Rubber Stand



A23.1503 Zoom Stereo Microscope, 0.7-4.5x, Mobile Phone Repair



A23.1503 With A59.4231 5.0M HDMI Camera



A23.1503 With A59.5107 7" LCD Camera



Real Time Display On Computer

Trinocular Head, Both Eye Tube Diopter Adjustable
Light Split Switch E100:P0/E0:P100
With Eyepiece Adapter + 1.0x C-Mount



Trinocular Head With CCD Adapter
Wide Field Eyepiece WF10x/22mm



Zoom Lens: 0.7~4.5x



Zoom Lens 0.7x~4.5x

Focusing Range 50mm, Distance To Center 155mm



Focusing Range 50mm, Distance To Center 155mm

Total Magnification Up To 3.5x-180x With Optional Eyepieces And Auxiliary Lens



**Total Magnification Up To 3.5x~180x
With Auxiliary Lens & Eyepiece**



Upper 56 LEDs Right Light Brightness Adjustable




Upper 56 LEDs Right Light Brightness Adjustable

Large Pole Stand, Size 480x320mm, Pole Height 240mm, Dia.32mm, With Heat Resist Soft Rubber Pad, Working Area For Phone PCB Repair





Large Working Area For Phone PCB Repair
Heat Resist Soft Rubber Pad





**Speical Designed Rubber Stand
For Mobile Phone PCB Repair Work**



**Speical Designed Rubber Stand
For Connecting Different Microscope**

A23.1503 DIY Phone PCB Repair Microscope	
Head	Trinocular Head, Both Eye Tube Diopter Adjustable Light Split Switch E100:P0/E0:P100 With Eyepiece Adapter + 1.0x C-Mount
Eyepiece	WF10x/20mm, Both With Eye-cup
Zoom Lens	0.7~4.5x
Zoom Ratio	1:6.5

Magnification	7x-45x, Up To 3.5x-180x With Optional Eyepieces And Auxiliary Lens	
W.D.	Working Distance 100mm, With Standard 1.0x Lens	
F.O.V	Field Of View 1~95.2mm	
Stand	Large Pole Stand, Size 480x320mm , Pole Height 240mm, Dia.32mm, With Heat Resist Soft Rubber Pad, Working Area For Phone PCB Repair	
Optional Accessories		Cata.No.
Eyepiece	WF15x/15mm, High-eyepoint	A51.1520-15
	WF20x/10mm, High-eyepoint	A51.1520-20
Auxiliary Lens	0.5x,W.D. 177mm	A52.1520-05
	0.75x,W.D. 120mm	A52.1520-75
	2.0x,W.D. 30mm	A52.1520-20
Light Source	Upper 56 LEDs Right Light Brightness Adjustable	A56.2103-A
Digital Camera	HDMI+USB+TF Digital Camera, 38.0M, USB Measure	A59.4231
	7" LCD Digital Camera, AV Output	A59.5107
	9" LCD Digital Camera, AV Output	A59.5109
CCD Adapter	Eyepiece Adapter + 0.5x C-Mount	A55.1520-05

Product Accessories



A51.1520-10
WF10x/20mm



A51.1520-15
WF15x/15mm



A51.1520-20
WF20x/10mm



A52.1520-05
0.5x,W.D.177mm



A52.1520-75
0.75x,W.D. 120mm



A52.1520-20
2.0x,W.D. 30mm



A56.2103-A Right Light



A59.4231 HDMI Camera



A55.1520-05
CCD Adapter



A59.5107 7" LCD



A59.5109 9" LCD

A23.1503 Zoom Stereo Microscope Optical Data Sheet								
Objective	Standard 1.0x		Auxiliary 0.5x		Auxiliary 0.75x		Auxiliary 2.0x	
W.D.	100mm		177mm		120mm		30mm	
Eyepiece	Mag	View Field	Mag	View Field	Mag	View Field	Mag	View Field
10X/20mm	7.0X	28.6mm	3.5X	57.1mm	5.25x	43.8mm	14X	14.3mm
	45.0X	4.4mm	22.5X	8.9mm	45x	6.5mm	90X	2.2mm
15X/15mm	10.5X	21.1mm	5.25X	42.8mm	7.87x	31.8mm	21X	10.7mm
	67.5X	3.3mm	33.75X	6.7mm	50.6x	4.7mm	135X	1.7mm
20X/10mm	14.0X	14.3mm	7.0X	28.6mm	10.5x	23.9mm	28X	7.1mm

	90.0X	2.2mm	45.0X	4.4mm	67.5x	3.6mm	180X	1.1mm
--	-------	-------	-------	-------	-------	-------	------	-------



Opto-Edu (Beijing) Co., Ltd.



0086 13911110627



sale@optoedu.com



cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China