



A12.0806 OPTO EDU Heating Stage Biological Microscope 195 Plan Achromatic Objective

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC
- Certification: CE
- Model Number: A12.0806
- Minimum Order Quantity: 1
- Price: \$1-\$1000
- Packaging Details: Carton
- Delivery Time: 1 month
- Payment Terms: T/T, West Union and Paypal
- Supply Ability: 10000



Product Specification

- Eyepiece: WF10x/20mm
- Nosepiece: 4 Holes Revolving Nosepiece, Inward
- Objective: 195 Plan Achromatic Objective
- Working Stage: Double Layer Mechanical Stage, Size: 160x142mm
- Highlight: **195 plan achromatic objective biological microscope**
, wf10x eyepiece biological microscope,
opto edu microscope for laboratory use



Product Description

Binocular & Trinocular, WF10x/20mm, WF16x/15, Quaduple, Plan 4x10x40x100x, Mechanical Stage 160x142/76x52mm, With Heater Temperature 20-50°, Heating Electric Power 12V36W, Abbe N.A.1.25, 3W LED



OPTO-EDU (BEIJING) CO., LTD.

OPTO-EDU

F-1501 Wanda Plaza, No.18 Shijingshan Road, Beijing 100043, China
Tel: +8610 88696020 Fax: +8610 88696085

A12.0806 Heating Stage Biological Microscope



Product Details



A12.0806

Heating Stage Biological Microscope

Heating Mechanical Stage

With Heater Temperature 20-50° ,
Heating Electric Power 12V36W,
Abbe N.A.1.25, 3W LED.

Easy To Operate

Built-in Handle For Easy Pick Up.

Product Outlook



Product Specification

A12.0806 Heating Stage Biological Microscope Specification		PB	PT
Head	Compensation Binocular Head, 30° Inclined, Interpupillary Distance 48-75mm	●	
	Compensation Trinocular Head, 30° Inclined, Interpupillary Distance 48-75mm		●
Eyepiece	WF10X/20mm	●	●
	WF16X/15mm	●	●
195 Plan Achromatic Objective	Plan 4x	●	●
	Plan 10x	●	●
	Plan 40x (s)	●	●

	Plan 100x (S, Oil)	●	●
Nosepiece	4 Holes Revolving Nosepiece, Inward	●	●
Focusing System	Coaxial Coarse & Fine Focusing Knob, Adjustable, Coarse Focusing Range 22mm, Fine Focusing Graduation 0.002mm	●	●
Stage	Double Layer Mechanical Stage, Size: 160x142mm Moving Range: 76x52mm Can Hold 2 Slides Temperature Control Range: Room Temperature 20°C-50°C Temperature Control Accuracy: $\pm 1^{\circ}\text{C}$ Heating Power: 12V Power: 36W	●	●
Condenser	Abbe Condenser N.A.=1.25, Iris Diaphragm, Center Adjustable, With Height Adjust Knob	●	●
Power	Switch Power, Wide Voltage AC90-240V	●	●
Light Source	Bottom Light, LED 3W/1A, Stepless Adjustable Light, No Screen Flashing	●	●
Dimming System	Digital Circuits, Precision Stepless Adjustable Light, No Screen Flash	●	●
Packing	Foam + Carton	●	●
Adapter	0.5x C-Mount	○	○
	0.65x C-Mount	○	○
	1x C-Mount	○	○
Note: "●" In Table Is Standard outfits, "○" Is Optional Accessories.			



Opto-Edu (Beijing) Co., Ltd.



0086 13911110627



sale@optoedu.com



cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China