



0.28-1.875x Stereo Optical Microscope 3d Motorized Rotate Led Illumination Zoom Mono A21.1609

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE,
- Model Number: A21.1609
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~3000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Zoom Range: 0.28-1.875x
- Illuminators: Long Life LED Illumination
- Miniature Motor: 90V-260V Wide Input Voltage, 5-27turns/minute, Stepless Adjustment
- Holder And Main Body Matching Dimension: 45mm
- CCD Adapter: 0.5x
- 1/3"22CCD FOV(mm): 12.9x17.1-1.92x2.56
- Highlight: **1.875x stereo optical microscope, stereo optical microscope 3d, 0.28x microscope with video output**

Product Description

Advanced Optical System & Precision Machinery Electronics Technology

3D Whole Azimuth Multi-Angle Observation Greatly Improved Detail Observation Level

360° Observe Each Side of Object With Large Depth of Field In Real Time Dynamic Image

Observe The Position Which Cannot Be Viewed By Traditional Lens Like Bottom, Inside, Side-Wall Of Holes


Widely Used For SMT, PCB, BGA Surface Board Assembly Industry



A21.1609 3D Rotated Zoom Video Microscope


45°/60° Optional, Lens 0.28x-1.875x, W.D. 26mm/45°, 42mm/60°, LED Light, Mini Motor 5-27 Turns/Min, Stepless Adjustment, Stand SD11-D, Moving Stage Range 100x86mm, 0.5x CCD Adapter, Control Box





A21.1609 3D Rotated Zoom Video Microscope		
Zoom Range	0.28-1.875x	
Observation Angle	45°	60°
Standard WD	26mm	42mm
Illuminators	Long Life LED Illumination	
Miniature Motor	90V-260V Wide Input Voltage, 5-27turns/minute, Stepless Adjustment	
Holder And Main Body Matching Dimension	45mm	
CCD Adapter	0.5x	
1/3"CCD FOV(mm)	12.9x17.1-1.92x2.56	

Specification



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China