



30 Fps Opto Edu A33.5100 7 Inch Digital Microscope Heating Stage Biological

Our Product Introduction

for more products please visit us on cnoec.com

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A33.5100
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month

OPTO-EDU



Product Specification

- Frame Rate: 30 Fps
- Output: USB 2.0, Output To LCD, PC Or Mobile Phone
- Power: DC5V/2A, Power Port 3.5/1.35 DC
- Storage: Micro SD Card Slot, Support Class 10 Max 32G
- Head: Teaching Head, 45°+90° Inclined, 360° Rotatable
- Eyepiece: WF10x, WF16x
- Highlight: opto edu 7 inch digital microscope, 30 fps 7 inch digital microscope, opto edu portable lcd digital microscope

OPTO-EDU




More Images


OPTO-EDU



OPTO-EDU



LED Illumination



Heated Stage



Heating Set

Product Description

7" LCD Monitor + Digital Camera All-in-One For Easy Use
Heating Stage Temperature Adjustable Range 30°~50° For Blood, Sperm View
Total Optical Magnification 40x~640x, Digital Magnification 200x~2000x On Screen
Coaxial Coarse & Fine Focusing Knob, Focus Range 20mm
Transmit Light Source LED Brightness Powered By 12V DC Adapter



OPTO-EDU (BEIJING) CO., LTD.

F-1501 Wanda Plaza, No.18 Shijingshan Road, Beijing 100043, China
Tel: +8610 88696020, E-mail: sale@cnoec.com

OPTO-EDU

A33.5100
7" LCD Digital Heating Stage
Biological Microscope

for more products please visit us on cnoec.com



**7" LCD Monitor + Digital Camera
All-in-One For Easy Use**







A33.5100 7" LCD Heating Stage Biological Microscope	
LCD & Digital Camera	
LCD	1/2.75" CMOS
Sensor	1/2.75" CMOS
Resolution	2.0M 1920x1080p, With 0.3x C-Mount
Pixel Size	2.8x2.8um
Frame Rate	30 fps
Output	USB 2.0, Output to LCD, PC or Mobile Phone
Power	DC5V/2A, Power Port 3.5/1.35 DC
Storage	Micro SD Card Slot, Support Class 10 Max 32G
LCD Menu	One Key Switch Image From Main Camera & Hand-held Camera
	One Key Take Photo, Video, Record, Save to SD Card
	Set Color, Contrast, White Balance Auto/Manual Adjust
	Set Photo/ Video Resolution Adjust, Mic On/Off to Record Voice With Video
	Set Date, Time
PC Software	Choose Language: Chinese, Traditional Chinese, English, Japanese, France, Germany, Korea, Russia, Italy, Portuguese, Spanish, Arabic
	Playback Photo, Video
	CD With Measuring Software Strange View, Support Win/iMac Computer, Support Language Chinese, English, Russia, Ukrian
Microscope	
Head	Teaching Head, 45°+90° Inclined, 360° Rotatable
Eyeiece	WF10x, WF16x
Nosepiece	Triple Nosepiece, Ball Bearing Inner Locating
Objective	Achromatic 4x,10x,40x(s)
Working Stage	Heated Stage 37°C-38°C, Size 115mmx125mm, Moving 60mmx30mm
Focusing	Coaxial Coarse & Fine Focus Knob, Coarse Focusing Range 20mm
Condenser	Abbe Condenser With Iris Diaphragm
Light Source	LED, Brightness Adjustable, Powered By 5V DC Adapter
Heating Set	S Key:
	Long Press: On/Off
	Short Press: Set Temp.
Accessory	Prepared Slide 10pcs in Box,
	Experiment Kit,
	Mobile Phone Connector
	Power Adapter
	Power Cord
Packing	USB Cable
	Software CD
Aluminium Box + Master Carton	



Opto-Edu (Beijing) Co., Ltd.



0086 13911110627



sale@optoedu.com



cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China