



Portable Cordless Biological Microscope 1000X For School / Home A11.1123

Our Product Introduction

for more products please visit us on cnoec.com

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A11.1123
- Minimum Order Quantity: 1 pc
- Price: Negotiation
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month




Product Specification

- Applications: Student, School, Hobby, Home
- Certification: CE, Rohs
- Drawtube: Binocular
- Illumination: LED
- Max Magnification: 1000X
- Product Name: Microscope
- Theory: Biological Microscope
- Highlight: **high powered microscope, binocular biological microscope**



More Images



Product Description

Portable Cordless Biological Microscope 1000X For School / Home A11.1123

Description:

We have various Student Biological Microscope for your choice for student educational purpose. Designed with high quality optical system, steady and heavy materials, we can supply a real optical microscope for all students at very good prices. Sold to world market in large quantity, these student microscopes has become most popular tools in primary, middle and high school.

Specifications:


A 11.1123 Economical Binocular Biological microscope

This scope is an ideal school microscope, it is with LED lighting and fully portable with no need for messy classroom wires! Arm with handle hole is easy to carry

Our Product Introduction

A11.1123 Economical Binocular Biological Microscope Specification	
Head	Sliding Binocular Head, 45°Inclined, 360°Rotatable
Eyeiece	WF10X with pointer
Objective	Achromatic 4X
	Achromatic 10X
	Achromatic 40X, Spring
Nosepiece	Achromatic 100X, Spring/Oil
	Quadruple Revolving Nosepiece
Stage	Double layer Mechanical Stage ,Size 100x90mm
Condenser	NA1.25 Condensor with Iris Diaphragm
Focusing	Coaxial Coarse and Fine Focusing Knobs
Light Source	Top and Bottom LED light with rechargeable function (if shipped by air, battery will not be included)
	Power supply adapter







Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China