



100X Hobby High Power Monocular Biological Microscope LED Illumination A11.1124

Our Product Introduction

for more products please visit us on cnoec.com

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A11.1124
- Minimum Order Quantity: 1 pc
- Price: Negotiation
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Applications: Student, School, Hobby, Home
- Certification: CE, Rohs
- Drawtube: Monocular
- Illumination: LED
- Product Name: Microscope
- Theory: Biological Microscope
- Highlight: binocular biological microscope, biological compound microscope



More Images



Product Description

High Power Monocular Biological Microscope LED Illumination A11.1124

Description:

- 40x~400x Magnification
- Up & Bottom LED Light Source
- Full Metal Body for Steady Observation
- Compact All in One Design for Easy Moving

Specifications:

Our Product Introduction

A11.1124 Economical Biological Microscope Specification	
Head	Monocular Head
Eyeiece	WF10X
Objective	Achromatic 4X
	Achromatic 10X
Nosepiece	Achromatic 40X, Spring
	Triple Revolving Nosepiece
Stage	Plain Stage with Clips
Condenser	Disc Diaphragm
Focusing	Coarse Focusing Knobs
Illumination	Top and Bottom LED light, powered by AA 1.5V Battery 3 pcs. (No battery inclued when ship by air)

About us :

OPTO-EDU have been focusing in microscope and educational instrument export market for many years. Our customers are from U.S.A, U.K., Russia, Canada, Germany, Denmark, Poland, Sweden, U.A.E. Qatar, Saudia Arabia, Egypt, Argentina, Korea, Thailand, Brazil, and more countries. Concentrating on Optical and Educational products, we are committed to be clients' excellent consultant and supplier for products and services in this field. We strive to be a unique, long-term and excellent supplier for all kinds of microscope from China. We sincerely welcome customers from all over the world contact us at any time!







Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China