



WF10x 100X Elementary Biological Microscope Research , Rohs CE A11.2601

Our Product Introduction

for more products please visit us on cnoec.com

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A11.2601
- Minimum Order Quantity: 1 pc
- Price: Negotiation
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Applications: Student, School
- Certification: CE, Rohs
- Drawtube: Binocular, Monocular
- Illumination: Halogen Lamp, LED
- Max Magnification: 1000X
- Product Name: Microscope
- Theory: Biological Microscope
- Highlight: **high powered microscope ,
biological compound microscope**



More Images



Product Description

1000X Elementary Biological Microscope Research , Rohs CE A11.2601

Description:

A11.2601 Student Microscope, modern design with high quality optical system, supplying clear image, is good to be used in high school, university for student.

- Anti-fugal optical system
- Wide voltage 80-265V, for all countries
- Built-in blue filter
- LED light source available

Specifications:

Specification	A11.2601-A	A11.2601-B	A11.2601-C
Head	Monocular head, inclined 30°, rotatable 360°	Binocular head, inclined 30°, rotatable 360° Interpupillary distance 50-75mm, Diopter adjustable ±5	
Eyepiece	WF10x/18mm		
Objective	Achromatic 4x/0.1, 10x/0.25, 40x/0.65(S), 100x/1.25(O,S)		Plan 4x/0.1, 10x/0.25, 40x/0.65(S), 100x/1.25(O,S)
Nosepiece	Quadruple nosepiece		
Focusing	Coaxial coarse and fine focusing knob, Coarse stroke 22mm, Fine division 2µm, Condenser up-down range 25mm		
Stage	Double layer mechanical, size 120X115mm, moving range 76 x50mm		
Condenser	Abbe N.A.1.25 with aperture iris diaphragm, Aperture center is adjustable		
Filter	Built-in blue filter Φ38.5mm		
Light source	6V 20W halogen lamp with brightness control(100V-240V)		
Packing	1pc/carton, size 34x24x42cm, G/N weight 6.5/5.5 Kgs		
Optional Accessories			Item No.
Head	Upgrade From Monocular to Binocular		A53.2651-B
	Upgrade From Binocular to Trinocular		A53.2652-B
Eyepiece	WF5x/20mm		A51.2602-05
	WF16x/13mm		A51.2602-16
	WF10x/18mm with pointer		A51.2604-10
	WF10x/18mm with micrometer 0.1mm/Div		A51.2605-10
	WF16x/13mm with micrometer 0.1mm/Div		A51.2605-16
Objective	Achromatic 20x		A52.2601-B20
	Achromatic 60x(S)		A52.2601-B60
	Plan A chromatic 20x		A52.2603-B20
	Plan Achromatic 60x(S)		A52.2603-B60
Light Source	Plano-concave mirror Φ52mm		A56.2611
	Upgrade to 3W LED, non-rechargeable		A56.2651-3W
	Upgrade to 0.3W LED, rechargeable		A56.2651-0.3W
Dark Field	Dark field ring for for 4x- 40x objective		A5D.2601-B
Adapter	CCD Adapter 0.5x (Focus Adjustable)		A55.2601-B05
	CCD Adapter 1.0x (Focus Adjustable)		A55.2601-B10
Polarizing Kit	Analyzer & Polarizer		A5P.2601-B





Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China