



10X Clinic Stereo Optical Microscope Low Magnification Microscopes A22.1204

Our Product Introduction

for more products please visit us on cnoec.com

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A22.1204
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Magnification: 10X
- Applications: Medical, Clinic
- Certification: CE, Rohs
- Drawtube: Binocular
- Product Name: Stereo Zoom Microscope
- Theory: Stereo Microscope
- Highlight: **stereo dissecting microscope, stereo binocular microscopes**



More Images



Product Description

20X Clinic Stereo Optical Microscope Low Magnification Microscopes A22.1204

Description:

Stereo Optical Microscope, designed for low magnification observation of a sample, generate three-dimensional, laterally precise images, to see where a three-dimensional view is needed for analyzing the detail. It is ideal for observing large biological samples or analysis of rough materials surfaces, widely used in industry production line, coin & jewelry inspection and etc.

Specifications:

A22.1204 Specification									
Model	Eyepiece	Objective	Working Stage	Illuminator					
1A	WF10X	Inserted 2x	White/Black	No					
1B	WF10X	Inserted 2x	White/Black	Incident					
1C	WF10X	Inserted 2x	Frosted Glass	Incident & Transmitted					
Optical Accessories				Item No.					
Eyepiece	WF5X			A51.1201-05					
	WF10X			A51.1201-10					
	WF15X			A51.1201-15					
	WF20X			A51.1201-20					
Objective	1x			A52.1280-01					
	2x			A52.1280-02					
	3x			A52.1280-03					
	4x			A52.1280-04					
	6x			A52.1280-06					
Technical Data									
Objective	Eyepiece								Working Distance (mm)
	Magnification	Field of View (mm)	Magnification	Field of View (mm)	Magnification	Field of View (mm)	Magnification	Field of View (mm)	
1X	5X	20.0	10X	20.0	15X	15.0	20X	10.0	57
2X	10X	10.0	20X	10.0	30X	7.5	40X	6.5	80
3X	15X	7.0	30X	6.7	45X	5.0	60X	4.3	61
4X	20X	5.0	40X	5.0	60X	3.7	80X	3.2	57
6X	30X	4.0	60X	3.3	90X	2.5	120X	2.1	66



A22.1204-1A

www.cnoec.com



A22.1204-1B



A22.1204-1C

www.cnoec.com





Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627 ✉ sale@optoedu.com 🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China