



Binocular Digital Stereo Microscope biological Microscopes 5x / 10x Rohs A22.1209

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A22.1209
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Magnification: 5-80X
- Applications: Medical, Clinic
- Certification: CE, Rohs
- Drawtube: Binocular
- Product Name: Stereo Zoom Microscope
- Theory: Stereo Microscope
- Highlight: stereo zoom microscopes,
stereo binocular microscopes



More Images



Product Description

Binocular Digital Stereo Microscope biological Microscopes Rohs A22.1209

Description:

Stereo Optical Microscope, designed for low magnification observation of a sample, generate three-dimensional, laterally precise images, to see where a three-dimensional view is needed for analyzing the detail. It is ideal for observing large biological samples or analysis of rough materials surfaces, widely used in industry production line, coin & jewelry inspection and etc.

Specifications:

A22.1209 Specification										
Model		Stand		Objective		Plate		Illuminator		
6A		Pole		Rotated 1x- 2x 1x-3x 2x-4x		White/Black		No		
6A-W		Arm Curved				White/Black		No		
6C		Pole				Frosted Glass		Incident & Transmitted		
6C-W		Arm Curved				Frosted Glass		Incident & Transmitted		
Technical Data										
Objective	Eyepiece								Working Distance (mm)	
	5X		10X		15X		20X			
	Magnification	Field of View (mm)	Magnification	Field of View (mm)	Magnification	Field of View (mm)	Magnification	Field of View (mm)		
	1X	5X	20.0	10X	20.0	15X	15.0	20X		10.0
	2X	10X	10.0	20X	10.0	30X	7.5	40X		6.5
	1X	5X	20.0	10X	20.0	15X	15.0	20X		10.0
	3X	15X	7.0	30X	6.7	45X	5.0	60X		4.3
2X	10X	10.0	20X	10.0	30X	7.5	40X	6.5		
4X	20X	5.0	40X	5.0	60X	3.7	80X	3.2		



A22.1209-6A

A22.1209-6AW



A22.1209-6C

A22.1209-6CW





Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China